

Apelin-36 (Rat, Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-057-28	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

[Apelin - Research Abstracts, Antibody Stainings, Kits Standard Curve, Sequences, etc.](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-057-28

Standard Size 96 wells


Sensitivity Minimum Detectable Concentration = 5.7 pg/ml

Range 0-10000 pg/ml

Linear Range 5.7 - 346 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Apelin-36 (Rat, Mouse)	100
	Apelin-12 (Human, Rat, Mouse Bovine)	100
	Apelin-13 (Human)	100
	Apelin-36 (Human)	100
	Alpha-ANP (Human)	0
	ANP (25-56) (N-ANP)	0
	ADM-52 (Human)	0
	BNP-32 (Human)	0
	CNP-22	0
	Ghrelin (Human)	0
	ET-1 (Human)	0
	Bradykinin	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$