

PYY (3-36) (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-059-02	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References

[PYY](#)

[PYY Immunoassay Kits](#)

[PYY Sequence Comparison](#)

[Normal Plasma Peptide/Protein Levels](#)

[PYY Antibodies](#)

[PYY EIA Kit References](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-059-02

Standard Size 96 wells


Sensitivity Minimum Detectable Concentration = 12 pg/ml

Range 0-10000 pg/ml

Linear Range 12-133 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Peptide YY (3-36) (Human)	100
	Peptide YY (Human)	100
	Peptide YY (Rat, Mouse, Porcine)	5
	Peptide YY (3-36)(Rat, Mouse, Porcine)	5
	Nuropeptide Y (Human)	0
	Neuropeptide Y (Porcine)	0
	Pancreatic Polypeptide (Human)	0
	Insulin (Human)	0
	Glucagon (Human)	0
	Amylin Amide (Human)	0
	Substance P	0

Standard
Curve Image



Please click [here](#) or click the reference tab on the top to see more articles that use this product.

Sample Type Human Extracted* Plasma (Pooled)

Linear Range ** 16.3-163 pg/ml

ED50 45 pg/ml

Sample Level **52.87 pg/ml**