

# PYY (Rat, Mouse, Porcine, Canine) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-059-03	1 kit	\$460.00	<input type="button" value="ORDER"/>

- [Product description](#)
- [More information](#)

- [PYY](#)
- [PYY Immunoassay Kits](#)
- [PYY Sequence Comparison](#)
- [Normal Plasma Peptide/Protein Levels](#)
- [PYY Antibodies](#)

**Protocol** [Fluorescent EIA Kit \(PDF\)](#)

**Sample Preparation** [Blood collection](#)  
[Plasma extraction](#)  
[Tissue extraction](#)  
[CSF extraction](#)

**Catalog #** FEK-059-03

**Standard Size** 96 well plate


**Sensitivity** Minimum Detectable Concentration = 16.2 pg/mL

**Range** 0-10000 pg/ml

**Linear Range** 16.2-227 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	Peptide YY (1-36) (Porcine, Rat)	100
	Peptide YY (3-36) (Porcine, Rat)	100
	Peptide YY (3-36) (Human)	3
	Peptide YY (1-36) (Human)	1

**Standard Curve Image**



4-P Fit:  $y = (A - D) / (1 + (x/C)^B) + D$

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.