

Vasopressin (AVP) [Arg8] (Human, Rat, Mouse, Ovine, Canine) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-065-07	1 kit	\$460.00	ORDER

Product description

More information


References

Normal Plasma Peptide/Protein Levels

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-065-07
Standard Size	96 well plate
Sensitivity	Minimum Detectable Concentration = 11.2 pg/mL
Range	0-10000 pg/ml
Linear Range	11.2-417 pg/ml

Cross Reactivity	Peptide	% Cross-reactivity
	[Arg ⁸]-Vasopressin (Human, Rat, Mouse, Bovine, Ovine, Canine)	100
	[Arg ⁸]-Vasotocin	78
	[Lys ⁸]-Vasopressin (Human, Rat, Mouse, Ovine)	0
	Deamino-[D-Arg ⁸]-Vasopressin	0
	Oxytocin	0
	LH-RH	0
	ACTH (Human)	0
	a-ANP(1-28) (Human, Canine)	0
	Met-Enkephalin	0
	PACAP38 (Human, Ovine, Rat)	0
	PACAP27-NH ₂ (Human, Ovine, Rat)	0

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$