

Secretin (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-067-03	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

References

Secretin

Normal Plasma Peptide/Protein Levels

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-067-03

Standard Size 96 well

Sensitivity Minimum Detectable Concentration = 57 pg/ml


Range 0-10000 pg/ml

Linear Range 57-1031 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Secretin (Human)	100
	Secretin (Porcine)	100
	Secretin (Rat)	100
	Secretin (Chicken)	0
	VIP (Human, Rat, Porcine)	0
	CCK (26-33) (Non Sulfated)	0
	Somatostatin	0
	NPY (Porcine)	0
	Met-Enkephalin	0
	Insulin (Human)	0
	Substance P	0
	PHI (Porcine)	0
	Amylin Amide (Human)	0

Sequence His - Ser - Asp - Gly - Thr - Phe - Thr - Ser - Glu - Leu - Ser - Arg - Leu - Arg - Glu - Gly - Ala - Arg - Leu - Gln - Arg - Leu - Leu - Gln - Gly - Leu - Val - NH₂

Standard
Curve Image



Please click [here](#) or click the reference tab on the top to see more articles that use this product.

Sample Type Human Extracted* Plasma (Pooled)

Linear Range ** 57-1031 pg/ml

ED50 168.61 pg/ml

Sample Level **13.5 pg/ml**