

CCK-33 (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-069-02	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

[All About Cholecystokinin](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-069-02

Standard Size 96 well

Sensitivity Minimum Detectable Concentration = 16.1 pg/ml

Range 0-10000 pg/ml

Linear Range 16.1-792 pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	CCK-33 (Human)	100
	CCK-33 (Human) (Non-Sulfated)	42.2
	CCK (26-33) (Non-Sulfated)	2.36
	CCK (1-21) (Human)	42.2
	CCK (1-21) (Porcine)	0
	CCK-8 (Non-Sulfated)	1
	CCK-8 (Human, Rat, Mouse) (26-33)	5.5
	Caerulein	0
	Gastrin-1 (Human)	0
	Pancreatic Polypeptide (Human)	0
	VIP (Human, Porcine, Rat)	0

Sequence Lys - Ala - Pro - Ser - Gly - Arg - Met - Ser - Ile - Val - Lys - Asn - Leu - Gln - Asn - Leu - Asp - Pro - Ser - His - Arg - Ile - Ser - Asp - Arg - Asp - Tyr(SO₃H) - Met - Gly - Trp - Met - Asp - Phe - NH₂

Standard
Curve Image

