

CCK (26-33) (Non-Sulfated) (Human, Rat, Mouse, Canine) - Fluorescent EIA Kit

| | | | |
|------------|---------------|----------|--------------------------------------|
| Catalog# | Standard size | Price | Order |
| FEK-069-04 | 1 kit | \$460.00 | <input type="button" value="ORDER"/> |

Product description

More information

References

[All About Cholecystokinin](#)

[Normal Plasma Peptide/Protein Levels](#)

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-069-04

Standard Size 96 well Fluorescent EIA Kit


Sensitivity Minimum Detectable Concentration = 15 pg/ml

Range 0-10000 pg/ml

Linear Range 15-333 pg/ml

| Cross Reactivity | Peptide | %Cross-reactivity |
|-------------------------|--|--------------------------|
| | CCK (26-33) (Non-Sulfated) (Human, Rat, Mouse) | 100 |
| | CCK-33 (Porcine) | 100 |
| | Caerulein | 100 |
| | Gastrin-1 (Human) | 100 |
| | Big Gastrin-1 (Human) | 100 |
| | CCK (27-33) | 100 |
| | CCK-33 (Non-Sulfated) (Human) | 77.8 |
| | CCK (Human) (1-21) | 75 |
| | CCK-33 (Sulfated) (Human) | 42.9 |
| | CCK (30-33) | 12.8 |
| | CCK (26-33) (Sulfated) | 3.9 |
| | Pancreatic Polypeptide (Human) | 0 |
| | VIP (Human, Porcine, Rat) | 0 |

Standard
Curve Image



4-P Fit: $y = (A - D) / (1 + (x/C)^B) + D$