


Urotensin II (Human) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-071-05	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description | **More information** | **References**

Urotensin II

Protocol	Fluorescent EIA Kit (PDF)
Sample Preparation	Blood collection Plasma extraction Tissue extraction CSF extraction
Catalog #	FEK-071-05
Standard Size	96 well
Sensitivity	Minimum Detectable Concentration = 8.4 pg/ml
Range	0-10000 pg/ml
Linear Range	8.4-843 pg/ml
Sequence	Glu - Thr - Pro - Asp - Cys - Phe - Trp - Lys - Tyr - Cys - Val ~ Disulfide Bridge between Cys5 & Cys10
Standard Curve Image	



$$4\text{-P Fit: } y = (A - D) / (1 + (x/C)^B) + D$$

Please click [here](#) or click the reference tab on the top to see more articles that use this product.

Sample Type	Human Extracted* Plasma (Pooled)
Linear Range **	8.4-843 pg/ml
ED50	70.08 pg/ml
Sample Level	9.3 pg/ml

Thursday 20 December, 2012

11984163 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.