

Urotensin II (Mouse) - Fluorescent EIA Kit

Catalog#	Standard size	Price	Order
FEK-071-08	1 kit	\$460.00	<input type="button" value="ORDER"/>

Product description

More information

Urotensin II

Protocol [Fluorescent EIA Kit \(PDF\)](#)

Sample Preparation [Blood collection](#)
[Plasma extraction](#)
[Tissue extraction](#)
[CSF extraction](#)

Catalog # FEK-071-08

Standard Size 96 well

Sensitivity Minimum Detectable Concentration = 7.9 pg/ml

Range 0-10000 pg/ml

Linear Range 7.9-249pg/ml

Cross Reactivity	Peptide	%Cross-reactivity
	Urotensin II (Mouse)	100
	Urotensin II (Human)	100
	Urotensin II (Rat)	100
	Endothelin-1	0
	Angiotensin II	0
	ADM (Rat)	0
	PAMP-20 (Human)	0
	alpha-ANP-28 (Human)	0
	CGRP (Rat)	0
	CNP-22	0
	BNP-25 (Rat)	0
	BNP-32 (Rat)	0
	Bradykinin	0

Sequence pGlu - His - Gly - Ala - Ala - Pro - Glu - Cys - Phe - Trp - Lys - Tyr - Cys - Ile

Disulfide Bridge: Cys8 - Cys13