

Obestatin (Rat, Mouse) - Purified IgG Antibody

Catalog#	Standard size	Price	Order
G-031-90	200 µg	\$495.00	<div>ORDER</div>
Product description	More information	References	

Obestatin - Physiological Opponent of Ghrelin

Catalog #	G-031-90
Standard Size	200 µg
Sequence	Phe - Asn - Ala - Pro - Phe - Asp - Val - Gly - Ile - Lys - Leu - Ser - Gly - Ala - Gln - Tyr - Gln - Gln - His - Gly - Arg - Ala - Leu - NH2
Species	Mouse, Rat
Host	Rabbit
Reconstitution	For consistent and reproducible results, reconstitute with 200µl of distilled water for the equivalent of undiluted antibody, immediately before use.
Storage Condition	For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.
Contents	Each vial contains 200µg of rabbit IgG lyophilized in 200µl PBS (137 mM NaCL, 2.7 mM KCL, 4.3 mM Na2HPO4, 1.4 mM KH2PO4, pH7.5).