

Neuronostatin-13 (Human, Porcine) - Purified IgG Antibody

Catalog#	Standard size	Price	Order
G-060-50	200 ug	\$450.00	<input type="button" value="ORDER"/>

Product description

More information

Neuronostatin

Catalog #	G-060-50
Family	Neuronostatin
Standard Size	200 ug
Species	Human, Porcine
Host	rabbit
Immunogen	human/porcine Neuronostatin-13
Storage Condition	For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.

Thursday 20 December, 2012

11986164 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.