

Neuronostatin-13 (Rat, Mouse) - Antibody for Immunohistochemistry

Catalog#	Standard size	Price	Order
H-060-48	100 ul	\$450.00	<input type="button" value="ORDER"/>

[Product description](#) [More information](#) [References](#)

Neuronostatin

Catalog #	H-060-48
Family	Neuronostatin
Standard Size	100 ul
Species	Rat, Mouse
Host	Rabbit
Immunogen	rat/mouse Neuronostatin-13
Storage Condition	For optimal results, use the antibody as soon as possible after reconstitution. Store in lyophilized form unless needed and reconstitute immediately before use. Once reconstituted, the antibody should be stable for a few days at -4°C. For storage up to a few months, prepare small aliquots after reconstitution and freeze at -20°C or -80°C. Repeated freeze thaw cycles should be strictly avoided.
Recommended Dilution Factor	1:200

Thursday 20 December, 2012

11986701 requests since Wednesday 05 April, 2006

Copyright © 2012 PHOENIX PHARMACEUTICALS, INC.