

Human SOX7 Antibody

Antigen Affinity-purified Polyclonal Goat IgG Catalog Number: AF2766

DESCRIPTION	
Species Reactivity	Human
Specificity	Detects human SOX7 in direct ELISAs and Western blots. In direct ELISAs and Western blots, less than 2% cross-reactivity with recombinant human SOX1, -2, -3, -9, -10, and -17 is observed.
Source	Polyclonal Goat IgG
Purification	Antigen Affinity-purified
Immunogen	E. coli-derived recombinant human SOX7 Ser249-Ala345 Accession # Q9BT81
Formulation	Lyophilized from a 0.2 µm filtered solution in PBS with Trehalose. See Certificate of Analysis for details. *Small pack size (-SP) is supplied as a 0.2 µm filtered solution in PBS.
APPLICATIONS	
Please Note: Optimal dilutions should be determined by each laboratory for each application. General Protocols are available in the Technical Information section on our website.	
	Recommended Sample Concentration
Western Blot	0.1 μg/mL Recombinant Human SOX7
Immunohistochemis	try 5-15 µg/mL Immersion fixed paraffin-embedded sections of human breast and lung
PREPARATION AND STORAGE	
Reconstitution	Reconstitute at 0.2 mg/mL in sterile PBS.
Shipping	The product is shipped at ambient temperature. Upon receipt, store it immediately at the temperature recommended below. *Small pack size (-SP) is shipped with polar packs. Upon receipt, store it immediately at -20 to -70 °C
Stability & Storage	Use a manual defrost freezer and avoid repeated freeze-thaw cycles. 12 months from date of receipt, -20 to -70 °C as supplied. 1 month, 2 to 8 °C under sterile conditions after reconstitution. 6 months, -20 to -70 °C under sterile conditions after reconstitution.

BACKGROUND

SOX7 belongs to the Sox (SRY-like HMG-box gene) family of transcription factors with diverse roles in development. SOX7 has been suggested to play a crucial role in parietal endoderm differentiation. SOX7 is expressed in many embryonic and adult tissues.

