

HDAC5(Ab-498) Antibody

Catalog No: #21142



Package Size: #21142-1 50ul #21142-2 100ul #21142-4 25ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | HDAC5(Ab-498) Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Purification | Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide. |
| Applications | WB IHC |
| Species Reactivity | Human Mouse Rat |
| Specificity | The antibody detects endogenous level of total HDAC5 protein. |
| Immunogen Type | Peptide-KLH |
| Immunogen Description | Peptide sequence around aa.496~500 (T-Q-S-S-P) derived from Human HDAC5/7. |
| Target Name | HDAC5 |
| Other Names | HD5 |
| Accession No. | Swiss-Prot: Q9UQL6NCBI Protein: NP_001015053.1 |
| Concentration | 1.0mg/ml |
| Formulation | Supplied at 1.0mg/mL in phosphate buffered saline (without Mg2+ and Ca2+), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol. |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |


Application Details

Predicted MW: 124kd
Western blotting: 1:500~1:1000
Immunohistochemistry: 1:50~1:100

Images



Western blot analysis of extracts from 293 and HepG2 cells using HDAC5(Ab-498) Antibody #21142.



Immunohistochemical analysis of paraffin-embedded human breast carcinoma tissue using HDAC5(Ab-498) Antibody #21142(left) or the same antibody preincubated with blocking peptide(right).

Background

Histones play a critical role in transcriptional regulation, cell cycle progression, and developmental events. Histone acetylation/deacetylation alters chromosome structure and affects transcription factor access to DNA. The protein encoded by HDAC5 belongs to the class II histone deacetylase/acuc/alpha family. It possesses histone deacetylase activity and represses transcription when tethered to a promoter. It coimmunoprecipitates only with HDAC3 family member and might form multicomplex proteins. It also interacts with myocyte enhancer factor-2 (MEF2) proteins, resulting in repression of MEF2-dependent genes. This gene is thought to be associated with colon cancer. Two transcript variants encoding different isoforms have been found for this gene.

Doppler H, et al. (2005) J Biol Chem. 280(15):15013-15019.

McKinsey TA, et al. (2000) Nature. 408(6808): 106-111.

Published Papers

Pang J, Yan C, Natarajan K et al., GIT1 mediates HDAC5 activation by angiotensin II in vascular smooth muscle cells., Arteriosclerosis, Thrombosis, and Vascular Biology, 28(5):892-898.(2008)

[PMID:18292392](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.