

CD45 Monoclonal Antibody

Catalog No: #40420



Package Size: #40420-1 50ul #40420-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD45 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous CD45 proteins.
Target Name	CD45
Other Names	PTPRC; CD45; Receptor-type tyrosine-protein phosphatase C; Leukocyte common antigen; L-CA; T200; CD45
Accession No.	Swiss-Prot: P08575NCBI Gene ID: 5788
SDS-PAGE MW	147kd
Concentration	1mg/ml
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:2000

Images



Western blot analysis of Hela, diluted at 1:2000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.