



Certificate of Analysis

Product: **AmpliStain™ anti-Goat 1-Step HRP**

Code #: **#AS-G1-HRP**

Description: **#AS-G1-HRP is a one step single reagent for the IHC staining and CISH with goat primary antibodies. #AS-G1-HRP comprises an enhanced anti-Goat IgG StrongZyme™ HRP conjugate diluted in Universal IHC Diluent/Blocker/Stabilizer (#UDBS) to the reserved final working strength**

Concentration: **Product supplied in Ready to Use RESERVE™ format. Further diluting with UDBS to the lower working strength concentrations may be recommended when detecting well accessible and/or abundant antigens**

Buffer base: **Universal IHC Diluent/Blocker/Stabilizer, #UDBS**

Lot: **192692**

Storage: **+2°C/+8°C. Do not freeze! / Do not store below 0°C!**

Expiration date: **05/2018**

Quality Control Assurance Release Procedure: **ELISA use test**

Result:

**POSITIVE
(COMPLIES)**

~~**NEGATIVE
(DOES NOT COMPLY)**~~

Date:

05.04.2017

This document has been electronically produced and is valid without a signature

< Strong, stronger, StrongZyme™ >>>