

ANTIBODIES FOR IMMUNOHISTO- CHEMISTRY

Catalog Number	Name	Size
H-006-23	ABC6 Protein (260-301) (Human) - Antibody for Immunohistochemistry	100 µl
H-072-13	Ac2-26 / Lipocortin-1 (2-26) / Annexin-1 (2-26) (Human) - Antibody for Immunohistochemistry	50 µl
H-001-06	ACTH (1-24) (Human, Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	50 µl
H-001-14	ACTH (18-39) / CLIP (Human) - Antibody for Immunohistochemistry	50 µl
H-001-01	ACTH (Human) - Antibody for Immunohistochemistry	50 µl
H-001-21	ACTH (Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-001-44	Adipo R1 (357-375) / Adiponectin Receptor Type 1 (357-375) (Human) - Antibody for Immunohistochemistry	50 µl
H-001-45	Adipo R1 (41-65) / Adiponectin Receptor Type 1 (41-65) (Human) - Antibody for Immunohistochemistry	50 µl
H-001-23	Adipo R2 (374-386) / Adiponectin Receptor Type 2 (374-386) (Human) - Antibody for Immunohistochemistry	50 µl
H-006-12	Adipocyte Fatty Acid Binding Protein / A-FABP (110-132) (Human) - Antibody for Immunohistochemistry	100 µl
H-ADI-01	Adiponectin Globular Form, Recombinant (Human) - Antibody for Immunohistochemistry	50 µl
H-ADI-06B	Adiponectin Globular Form, Recombinant (Mouse) - Antibody for Immunohistochemistry	500 µl
H-ADI-06A	Adiponectin Globular Form, Recombinant (Mouse) - Antibody for Immunohistochemistry	100 µl
H-ADI-06	Adiponectin Globular Form, Recombinant (Mouse) - Antibody for Immunohistochemistry	50 µl
H-028-52	Adiponectin Precursor (78-98) Amide (Human) - Antibody for Immunohistochemistry	200 µl
H-028-51	Adiponectin Precursor [Cys0] (222-244) C-terminal (Human) - Antibody for Immunohistochemistry	50 µl
H-028-54B	Adiponectin Precursor [Cys0] (223-244) (Rat, Mouse) - Antibody for Immunohistochemistry	200 µl
H-028-54A	Adiponectin Precursor [Cys0] (223-244) (Rat, Mouse) - Antibody for Immunohistochemistry	100 µl
H-002-72	Adiponutrin (378-413) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-028-47	Adiponutrin (436-481) (Human) - Antibody for Immunohistochemistry	100 µl
H-010-28	Adrenomedullin (1-44) (Human) - Antibody for Immunohistochemistry	100 µl
H-010-31	Adrenomedullin (AM / ADM) (1-50) (Mouse) - Antibody for Immunohistochemistry	50 µl
H-010-08	Adrenomedullin (AM / ADM) (1-50) (Rat) - Antibody for Immunohistochemistry	50 µl
H-010-01	Adrenomedullin (AM / ADM) (1-52) (Human) - Antibody for Immunohistochemistry	50 µl
H-010-07	Adrenomedullin (AM / ADM) (153-185), Pro / Adrenotensin (Human) - Antibody for Immunohistochemistry	50 µl
H-010-06	Adrenomedullin (AM / ADM) (45-92), Pro (Human) - Antibody for Immunohistochemistry	50 µl
H-010-24	Adrenomedullin N-20, Pro / PAMP-20 (Canine) - Antibody for Immunohistochemistry	50 µl
H-010-05	Adrenomedullin N-20, Pro / PAMP-20 / Prodepin (Human) - Antibody for Immunohistochemistry	50 µl
H-010-10	Adrenomedullin N-20, Pro / PAMP-20 / Prodepin (Rat) - Antibody for Immunohistochemistry	100 µl
H-010-40	Adrenomedullin-5 (Porcine) - Antibody for Immunohistochemistry	100 µl
H-010-37	Adrenomedullin-Gly (Human) - Antibody for Immunohistochemistry	100 µl
H-033-44	AF-16 Peptide (Human) - Antibody for Immunohistochemistry	100 µl
H-003-50	Agouti (1-40)-Amide (Human) - Antibody for Immunohistochemistry	50 µl

Catalog Number	Name	Size
H-003-55	AGRP (71-132) Form C-Amide (Human) - Antibody for Immunohistochemistry	50 µl
H-003-57	AGRP (82-131)-Amide (Mouse) - Antibody for Immunohistochemistry	50 µl
H-003-53	AGRP (83-132)-Amide (Human) - Antibody for Immunohistochemistry	50 µl
H-026-34	Alarin (Human) - Antibody for Immunohistochemistry	100 µl
H-026-33	Alarin (Rat) - Antibody for Immunohistochemistry	100 µl
H-017-06	Amylin (1-13) (Human) - Antibody for Immunohistochemistry	100 µl
H-017-01	Amylin (Feline) - Antibody for Immunohistochemistry	50 µl
H-017-03	Amylin (Human) - Antibody for Immunohistochemistry	50 µl
H-017-11	Amylin (Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-017-04	Amylin-Amide (Human) - Antibody for Immunohistochemistry	100 µl
H-018-01	Amyloid beta-Protein (1-40) (Human) - Antibody for Immunohistochemistry	50 µl
H-028-92	Angiopoietin-related Growth Factor (AGF) / Angiopoietin-like 6 (215-230) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-002-24	Angiotensin I (1-7) / Angiotensin II (1-7) (Human, Rat, Mouse, Canine) - Antibody for Immunohistochemistry	50 µl
H-002-01	Angiotensin I (Human, Rat, Mouse, Canine) - Antibody for Immunohistochemistry	50 µl
H-002-12	Angiotensin II (Human, Rat, Mouse, Canine) - Antibody for Immunohistochemistry	50 µl
H-005-21	ANP (104-116), Prepro (Human) - Antibody for Immunohistochemistry	50 µl
H-005-22	ANP (104-123), Prepro (Human) - Antibody for Immunohistochemistry	50 µl
H-005-20	ANP (56-92), Prepro (Human) - Antibody for Immunohistochemistry	50 µl
H-005-24	ANP (Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-005-06	ANP, alpha (1-28) (Human, Ovine, Canine) - Antibody for Immunohistochemistry	50 µl
PG-H-005-06	ANP, alpha (1-28) (Human, Ovine, Canine) - Antibody for Immunohistochemistry from Guinea Pig	50 µl
H-005-14	ANP, beta (Human) - Antibody for Immunohistochemistry	50 µl
H-001-79	Apelin Receptor APJ (357-377) (Rat) - Antibody for Immunohistochemistry	50 µl
H-057-23	Apelin-12 (Human, Rat, Mouse, Bovine) - Antibody for Immunohistochemistry	50 µl
H-057-15	Apelin-36 (Human) - Antibody for Immunohistochemistry	50 µl
H-057-28	Apelin-36 (Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-001-50	AUF-1 (151-199) (Human) - Antibody for Immunohistochemistry	50 µl
H-012-23	Augurin, Prepro (108-132) (Human) - Antibody for Immunohistochemistry	100 µl
H-012-24	Augurin, Prepro (133-148) (Human) - Antibody for Immunohistochemistry	100 µl
H-012-22	Augurin, Prepro (71-107) (Human) - Antibody for Immunohistochemistry	100 µl
H-008-29	Autotaxin (369-391) (Human, Mouse) - Antibody for Immunohistochemistry	100 µl
H-008-30	Autotaxin (754-783) (Human, Mouse) - Antibody for Immunohistochemistry	100 µl
H-048-61	AXOR12 (375-398) (Human) - Antibody for Immunohistochemistry	50 µl
H-030-51	BAFF-R [Cys0] (151-175) (Mouse) / BAFF-R [Cys0] (159-183) (Human) - Antibody for Immunohistochemistry	50 µl
H-072-50	Beacon (47-73) - Antibody for Immunohistochemistry	50 µl
H-011-26	BNP (1-45), Pro (Rat) - Antibody for Immunohistochemistry	50 µl
H-011-24	BNP (1-46), Pro (Human) - Antibody for Immunohistochemistry	50 µl
H-011-06	BNP (22-46), Pro (Human) - Antibody for Immunohistochemistry	50 µl
H-011-08	BNP-26 (Porcine, Ovine) - Antibody for Immunohistochemistry	50 µl
H-011-22	BNP-32 (Canine) - Antibody for Immunohistochemistry	50 µl
H-011-03	BNP-32 (Human) - Antibody for Immunohistochemistry	50 µl
H-011-10	BNP-32 (Porcine) - Antibody for Immunohistochemistry	50 µl
H-011-14	BNP-32 (Rat) - Antibody for Immunohistochemistry	50 µl
H-011-23	BNP-45 (Mouse) - Antibody for Immunohistochemistry	50 µl
H-011-17	BNP-45 (Rat) - Antibody for Immunohistochemistry	50 µl

Catalog Number	Name	Size
H-007-01	Bombesin - Antibody for Immunohistochemistry	50 µl
H-001-82	Bombesin Receptor Subtype-3 (BRS-3) (371-399) (Human) - Antibody for Immunohistochemistry	100 µl
H-009-38	Bombinakinin-GAP - Antibody for Immunohistochemistry	50 µl
H-009-01	Bradykinin (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-035-01	C-Peptide (Human) - Antibody for Immunohistochemistry	50 µl
H-035-03	C-Peptide 1 (Rat) - Antibody for Immunohistochemistry	100 µl
H-014-02	Calcitonin (Eel) - Antibody for Immunohistochemistry	50 µl
H-014-03	Calcitonin (Human) - Antibody for Immunohistochemistry	50 µl
H-014-06	Calcitonin (Rat) - Antibody for Immunohistochemistry	50 µl
H-014-09	Calcitonin (Salmon) - Antibody for Immunohistochemistry	50 µl
H-002-81	Cannabinoid Receptor CB2 (323-360) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-87	Carnitine Palmitoyl Transferase I, brain form (769-798) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-003-63	CART (1-39) (Human, Rat) - Antibody for Immunohistochemistry	50 µl
H-003-60	CART (55-102) (Human) - Antibody for Immunohistochemistry	50 µl
H-003-62	CART (55-102) (Rat, Mouse, Bovine) - Antibody for Immunohistochemistry	50 µl
H-003-61	CART (61-102) (Human, Rat, Mouse, Bovine) - Antibody for Immunohistochemistry	50 µl
H-044-01	Casomorphin beta - Antibody for Immunohistochemistry	50 µl
H-053-27	Catestatin (Human) - Antibody for Immunohistochemistry	100 µl
H-053-28	Catestatin (Mouse) - Antibody for Immunohistochemistry	100 µl
H-053-29	Catestatin (Rat) - Antibody for Immunohistochemistry	100 µl
H-002-82B	Cathepsin S (293-331) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-82	Cathepsin S (298-331) (Human) - Antibody for Immunohistochemistry	50 µl
H-002-63	Cathepsin S (34-57) (Human) - Antibody for Immunohistochemistry	100 µl
H-069-04	CCK (26-33) (Non-Sulfated) (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-069-27	CCK (46-69) (Human) - Antibody for Immunohistochemistry	100 µl
H-069-02	CCK-33 (Human) - Antibody for Immunohistochemistry	100 µl
H-015-02	CGRP (Human) - Antibody for Immunohistochemistry	50 µl
H-015-09	CGRP (Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-015-07	CGRP II (Human) - Antibody for Immunohistochemistry	50 µl
H-001-37	CGRP Type 1 Receptor (436-461) (Human) - Antibody for Immunohistochemistry	50 µl
H-002-51	Chemerin / TIG2 (135-154) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-50	Chemerin / TIG2 (145-157) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-52	Chemerin / TIG2 (29-67) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-53	Chemerin / TIG2 (72-93) (Human) - Antibody for Immunohistochemistry	100 µl
H-CHIP-01	CHIP (Human) - Antibody for Immunohistochemistry	100 µl
H-012-03	CNP-22 (Human, Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	50 µl
H-012-06	CNP-53 (1-29) (Rat, Porcine) - Antibody for Immunohistochemistry	50 µl
H-002-93	CNTF (27-72) (Human) - Antibody for Immunohistochemistry	100 µl
H-065-32B	Copeptin (Human) - Antibody for Immunohistochemistry	100 µl
H-065-32A	Copeptin (Human) - Antibody for Immunohistochemistry	50 µl
H-019-06	Corticotropin Releasing Factor (CRF) (Human, Rat, Mouse, Canine, Feline) - Antibody for Immunohistochemistry	50 µl
H-019-10	Corticotropin Releasing Factor (CRF) (Ovine) - Antibody for Immunohistochemistry	50 µl
H-006-24	Corticotropin Releasing Factor Receptor 2 (CRF-R2) [Cys0] (385-411) (Human) - Antibody for Immunohistochemistry	100 µl
H-060-04	Cortistatin-14 (Rat) - Antibody for Immunohistochemistry	50 µl
H-060-11	Cortistatin-17 (Human) - Antibody for Immunohistochemistry	50 µl
H-060-12	Cortistatin-29 (Human) - Antibody for Immunohistochemistry	50 µl
H-060-05	Cortistatin-29 (Rat) - Antibody for Immunohistochemistry	50 µl

Catalog Number	Name	Size
H-007-53	Coupling Factor 6 (CF6) Amide (Human) - Antibody for Immunohistochemistry	50 µl
H-007-54	Coupling Factor 6 (CF6) Precursor (55-108) Amide (Human) - Antibody for Immunohistochemistry	50 µl
H-072-15	CRAMP (140-173) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-015-64	CRSP-3 (Porcine) - Antibody for Immunohistochemistry	100 µl
H-072-53	Defensin 1, beta (Human) (synthetic) - Antibody for Immunohistochemistry	100 µl
H-072-48	Defensin 2, beta (Human) (synthetic) - Antibody for Immunohistochemistry	100 µl
H-072-42	Defensin 3, beta (Human) - Antibody for Immunohistochemistry	100 µl
H-072-44	Defensin 8, beta (Mouse) - Antibody for Immunohistochemistry	100 µl
H-028-82	Desnutrin (296-318) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-028-83	Desnutrin (465-486) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-013-01	DNP (Green Mamba) - Antibody for Immunohistochemistry	50 µl
H-021-21	Dynorphin A (1-13) (Human, Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	50 µl
H-021-10	Dynorphin A (1-8) (Human, Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	50 µl
H-021-03	Dynorphin A (Human, Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	50 µl
H-075-16	EC SOD (Extracellular Superoxide Dismutase) [Cys0] (33-54) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-047-69	EM-66 / Secretogranin II (219-284) (Rat) - Antibody for Immunohistochemistry	100 µl
H-044-10	Endomorphin-1 - Antibody for Immunohistochemistry	50 µl
H-044-11	Endomorphin-2 - Antibody for Immunohistochemistry	50 µl
H-021-40	Endorphin, alpha-Neo (Porcine) - Antibody for Immunohistochemistry	50 µl
H-022-06	Endorphin, beta (Bovine, Ovine, Camel) - Antibody for Immunohistochemistry	50 µl
H-022-14	Endorphin, beta (Human) - Antibody for Immunohistochemistry	50 µl
H-022-33	Endorphin, beta (Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-023-18	Endothelin-1 (109-125)-Amide, Prepro (Human) - Antibody for Immunohistochemistry	100 µl
H-023-05	Endothelin-1 (ET-1) (145-162)-Amide, Prepro (Human) - Antibody for Immunohistochemistry	100 µl
H-023-37	Endothelin-1 (ET-1) (18-50), Prepro (Human) - Antibody for Immunohistochemistry	50 µl
H-023-26	Endothelin-1 (ET-1) (94-109), Prepro (Porcine) - Antibody for Immunohistochemistry	50 µl
H-023-01	Endothelin-1 (ET-1) (Human, Rat, Mouse, Porcine, Bovine, Canine) - Antibody for Immunohistochemistry	50 µl
H-023-12	Endothelin-1, Big (22-38) (Human) - Antibody for Immunohistochemistry	50 µl
H-023-10	Endothelin-1, Big (Human) - Antibody for Immunohistochemistry	50 µl
H-023-13	Endothelin-2 (ET-2) (Human, Canine) - Antibody for Immunohistochemistry	50 µl
H-023-17	Endothelin-3 (ET-3) (Human, Rat, Porcine, Rabbit) - Antibody for Immunohistochemistry	50 µl
H-024-21	Enkephalin-Leucine - Antibody for Immunohistochemistry	50 µl
H-003-75	Enterostatin (APGPR) (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-070-94	Exendin-4 (Heloderma suspectum) - Antibody for Immunohistochemistry	100 µl
H-072-30	F2L / HBP (1-21) Nonacetylated (Human, Porcine, Canine) - Antibody for Immunohistochemistry	100 µl
H-072-34	F2L / HBP (1-21) Nonacetylated (Rat, Mouse) - Antibody for Immunohistochemistry	100 µl
H-001-65	Farnesoid X Receptor (FXR) (446-476) (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-028-77	FDC-SP (61-85) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-46	FGF-21 (26-47) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-48	FGF-21 (98-124) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-47	FGF-21 [Cys0] (184-209) (Human) - Antibody for Immunohistochemistry	100 µl
H-028-46	FIZZ-1 / RELM-beta (25-50) (Human) - Antibody for Immunohistochemistry	50 µl
H-028-44	FIZZ-2 / RELM-beta (25-45) (Mouse) - Antibody for Immunohistochemistry	50 µl
H-028-45	FIZZ-3 / Resistin (26-49) (Mouse) - Antibody for Immunohistochemistry	50 µl
H-047-29	FMRF-Amide (Molluscan Cardioexcitatory Neuropeptide) - Antibody for Immunohistochemistry	50 µl
H-026-01	Galanin (Human) - Antibody for Immunohistochemistry	50 µl
H-026-06	Galanin (Porcine) - Antibody for Immunohistochemistry	50 µl

Catalog Number	Name	Size
H-026-13	Galanin (Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-026-55	Galanin-Like Peptide (GALP) (36-60) (Human) - Antibody for Immunohistochemistry	50 µl
H-026-51	Galanin-Like Peptide (GALP) (Human) - Antibody for Immunohistochemistry	50 µl
H-026-52	Galanin-Like Peptide (GALP) (Rat) - Antibody for Immunohistochemistry	50 µl
H-027-04	Gastrin I (Human) - Antibody for Immunohistochemistry	50 µl
H-027-06	Gastrin I, Big (Human) - Antibody for Immunohistochemistry	50 µl
H-027-15	Gastrin Releasing Peptide (GRP) (1-16) (Porcine) - Antibody for Immunohistochemistry	50 µl
H-027-07	Gastrin Releasing Peptide (GRP) (Human) - Antibody for Immunohistochemistry	50 µl
H-027-13	Gastrin Releasing Peptide (GRP) (Porcine) - Antibody for Immunohistochemistry	50 µl
H-027-36	Gastrin Releasing Peptide (GRP), Pro (80-97) [Tyr0] (Human) - Antibody for Immunohistochemistry	100 µl
H-031-77	Ghrelin (52-75), Prepro (Human) - Antibody for Immunohistochemistry	50 µl
H-031-34	Ghrelin (86-117), Prepro (Human) - Antibody for Immunohistochemistry	50 µl
H-031-30	Ghrelin (Human) - Antibody for Immunohistochemistry	100 µl
H-031-52	Ghrelin (Porcine) - Antibody for Immunohistochemistry	100 µl
H-031-31	Ghrelin (Rat, Mouse) - Antibody for Immunohistochemistry	100 µl
H-032-11	Ghrelin O-acyltransferase (GOAT) (181-199) (Human)	100 µl
H-032-12	Ghrelin O-acyltransferase (GOAT) (356-375) (Human) - Antibody for Immunohistochemistry	100 µl
H-027-02	GIP (Human) - Antibody for Immunohistochemistry	50 µl
H-028-17	Glicentin-related Polypeptide (21-50) (Human) - Antibody for Immunohistochemistry	100 µl
H-028-26	GLP-1 (12-24) (Human) - Antibody for Immunohistochemistry	100 µl
H-028-25	GLP-1 (2-24) - Antibody for Immunohistochemistry	100 µl
H-028-11	GLP-1 (7-36)-Amide (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-028-13	GLP-1 (7-37) (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-028-14	GLP-2 [Arg34] (Human, Rat) - Antibody for Immunohistochemistry	50 µl
H-028-05	Glucagon (19-29) (Human, Rat, Mouse, Porcine, Bovine) - Antibody for Immunohistochemistry	100 µl
H-028-02	Glucagon (Human, Rat, Mouse, Porcine, Bovine) - Antibody for Immunohistochemistry	50 µl
H-029-07	Gn-RH Associated Peptide (GAP) (Rat) - Antibody for Immunohistochemistry	50 µl
H-002-83	GPR-30 (1-23) [Cys24] (Human) - Antibody for Immunohistochemistry	100 µl
H-002-84	GPR-30 (351-375) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-85	GPR-39 (423-453) (Human) - Antibody for Immunohistochemistry	100 µl
H-001-64	GPR-54 (375-396) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-031-02	Growth Hormone Releasing Factor (GHRF) (Human) - Antibody for Immunohistochemistry	50 µl
H-031-19	Growth Hormone Releasing Factor (GHRF) (Rat) - Antibody for Immunohistochemistry	50 µl
H-001-62	Growth Hormone Secretagogue Receptor type 1A [Cys0] (330-366) (Human) - Antibody for Immunohistochemistry	50 µl
H-001-61	Growth Hormone Secretagogue Receptor type 1B (245-289) (Human) - Antibody for Immunohistochemistry	50 µl
H-018-25	Humanin, Mitochondria-encoded / HN(M) (Human) - Antibody for Immunohistochemistry	50 µl
H-018-28	Humanin, Nuclear-encoded / HN(N) (Rat) - Antibody for Immunohistochemistry	50 µl
H-048-78	Hypothetical Protein XP_294524 (154-196) [pGlu154] Amide / QRFP-43 (Human) - Antibody for Immunohistochemistry	50 µl
H-035-10	IGF-1 (24-41) (Human) - Antibody for Immunohistochemistry	50 µl
H-035-15	IGF-1 (Human) - Antibody for Immunohistochemistry	50 µl
H-035-63	IGF-I, des-C-Peptide (Human) - Antibody for Immunohistochemistry	100 µl
H-034-01	Inhibin Alpha-Subunit (1-32) (Human) - Antibody for Immunohistochemistry	50 µl
H-034-05	Inhibin Alpha-Subunit (1-32) (Porcine) - Antibody for Immunohistochemistry	50 µl
H-035-43	INSL 3 (Mouse) - Antibody for Immunohistochemistry	100 µl
H-035-27	INSL 3 / RLF (Human) - Antibody for Immunohistochemistry	50 µl
H-035-48	INSL 3 C-Peptide (Human) - Antibody for Immunohistochemistry	100 µl
H-035-47	INSL 3 C-Peptide (Mouse) - Antibody for Immunohistochemistry	50 µl

Catalog Number	Name	Size
H-035-33	INSL 4 / RLF (Human) - Antibody for Immunohistochemistry	50 µl
H-035-56	INSL 4 Prepro-C-Peptide (64-110) (Human) - Antibody for Immunohistochemistry	100 µl
H-035-40	INSL 5 (Mouse) - Antibody for Immunohistochemistry	50 µl
H-035-58	INSL 5 C-Peptide, Prepro (49-106) (Human) - Antibody for Immunohistochemistry	100 µl
H-001-86	INSL 5 Receptor / GPCR-142 (310-330) (Human) - Antibody for Immunohistochemistry	100 µl
H-035-30	INSL 6 / RLF / RIF-1 (Human) - Antibody for Immunohistochemistry	50 µl
H-035-28	INSL 6 / RLF / RIF-1 A-Chain (Human) - Antibody for Immunohistochemistry	50 µl
H-035-29	INSL 6 / RLF / RIF-1 B-Chain (Human) - Antibody for Immunohistochemistry	50 µl
H-035-45	INSL 6 A-Chain (Mouse) - Antibody for Immunohistochemistry	50 µl
H-035-57	INSL 6 C-Peptide, Prepro (55-118) (Human) - Antibody for Immunohistochemistry	100 µl
H-035-55	INSL 6 Prepro-C-Peptide (121-148) (Human) - Antibody for Immunohistochemistry	100 µl
H-035-36	INSL 7 / Relaxin 3 (Human) - Antibody for Immunohistochemistry	50 µl
H-035-35	INSL 7 [Cys(Acm)10] / Relaxin 3 B-Chain [Cys(Acm)10] - Antibody for Immunohistochemistry	50 µl
H-035-52	INSL-7, C Peptide (Human) - Antibody for Immunohistochemistry	100 µl
PG-H-035-06	Insulin (Human) - Antibody for Immunohistochemistry	50 µl
H-035-53	Insulin-1 C-Peptide (Mouse) - Antibody for Immunohistochemistry	100 µl
H-028-85	Intelectin / Omentin (163-199) (Human) - Antibody for Immunohistochemistry	100 µl
H-028-87	Intelectin / Omentin (94-145) (Human) - Antibody for Immunohistochemistry	100 µl
H-010-57	Intermedin / Adrenomedullin-2 (17-47) (Human) - Antibody for Immunohistochemistry	50 µl
H-010-49	Intermedin / Adrenomedullin-2 (8-47) (Human) - Antibody for Immunohistochemistry	50 µl
H-010-50	Intermedin / Adrenomedullin-2 (Mouse) - Antibody for Immunohistochemistry	50 µl
H-010-52	Intermedin / Adrenomedullin-2 (Rat) - Antibody for Immunohistochemistry	50 µl
H-010-61	Intermedin / Adrenomedullin-2, Prepro (25-56) (Human) - Antibody for Immunohistochemistry	50 µl
H-010-62	Intermedin / Adrenomedullin-2, Prepro (57-92) (Human) - Antibody for Immunohistochemistry	50 µl
H-048-56	KISS-1 (112-121) Amide / Kisspeptin-10 / Metastin (45-54) Amide (Human) - Antibody for Immunohistochemistry	50 µl
H-048-59	KISS-1 (68-121) Amide / Metastin (1-54) Amide (Human) - Antibody for Immunohistochemistry	50 µl
H-048-62	KISS-1 (68-92) / Metastin (1-25) (Human) - Antibody for Immunohistochemistry	50 µl
H-048-60	KISS-1 (94-121)-Amide / Metastin (27-54)-Amide (Human) - Antibody for Immunohistochemistry	100 µl
H-005-19	LANP / ANP (26-55), prepro (Human) - Antibody for Immunohistochemistry	50 µl
H-075-40	LEAP-2 (38-77) (Human) - Antibody for Immunohistochemistry	100 µl
H-003-01	Leptin (116-167) (Human) - Antibody for Immunohistochemistry	50 µl
H-003-06	Leptin (22-56) (Human) - Antibody for Immunohistochemistry	50 µl
H-003-02	Leptin (57-92) (Human) - Antibody for Immunohistochemistry	50 µl
H-003-12	Leptin (Human) - Antibody for Immunohistochemistry	50 µl
H-003-13	Leptin (Mouse) - Antibody for Immunohistochemistry	100 µl
H-003-17	Leptin (Rat) - Antibody for Immunohistochemistry	100 µl
H-021-39	Leuomorphin / Enkephalin B (226-254), Prepro (Human) - Antibody for Immunohistochemistry	50 µl
H-001-53	LGR7 (735-757) (Human) - Antibody for Immunohistochemistry	50 µl
H-001-54	LGR8 (737-754) [Cys0] (Human) - Antibody for Immunohistochemistry	50 µl
H-040-02	LH-RH / Gn-RH (Human, Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	50 µl
H-075-06	LL37 (Human) - Antibody for Immunohistochemistry	50 µl
H-047-98	Manserlin / Sgll (529-568) (Rat) - Antibody for Immunohistochemistry	100 µl
H-006-17	MC4 Receptor (306-332) (Human) - Antibody for Immunohistochemistry	100 µl
H-006-16	MC4 Receptor (306-332) (Rat) - Antibody for Immunohistochemistry	100 µl
H-071-20	MCH-2R [Cys0] (237-251) (Human) - Antibody for Immunohistochemistry	50 µl
H-070-47	Melanin Concentrating Hormone (MCH) (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-033-35	MGF (Human) - Antibody for Immunohistochemistry	100 µl

Catalog Number	Name	Size
H-005-42	Mitocryptide-1 (MCT-1) (Human) - Antibody for Immunohistochemistry	100 µl
H-005-41	Mitocryptide-1 (MCT-1) (Porcine) - Antibody for Immunohistochemistry	100 µl
H-045-01	Motilin (Canine) - Antibody for Immunohistochemistry	50 µl
H-045-03	Motilin (Human, Porcine) - Antibody for Immunohistochemistry	50 µl
H-045-04	Motilin, Prepro (50-70) (Human) - Antibody for Immunohistochemistry	100 µl
H-045-06	Motilin, Prepro (93-115) (Human) - Antibody for Immunohistochemistry	100 µl
H-043-01	MSH, Alpha (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-043-12	MSH, Beta (Human) - Antibody for Immunohistochemistry	50 µl
H-043-16	MSH, Gamma (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-043-19	MSH, Gamma3 (Human) - Antibody for Immunohistochemistry	50 µl
H-028-65	Musclin / Osteocrin (Ostn) (113-130) (Mouse) / (115-132) (Rat) - Antibody for Immunohistochemistry	50 µl
H-028-72	Musclin / Osteocrin (Ostn) (31-77) (Rat) - Antibody for Immunohistochemistry	50 µl
H-028-70	Musclin / Osteocrin (Ostn) (31-78) (Human) - Antibody for Immunohistochemistry	50 µl
H-028-62	Musclin / Osteocrin (Ostn) (80-102) (Mouse) / (82-104) (Rat) - Antibody for Immunohistochemistry	50 µl
H-028-66	Musclin / Osteocrin (Ostn) (83-132)-Amide (Human) - Antibody for Immunohistochemistry	50 µl
H-028-68	Musclin / Osteocrin (Ostn) (83-133) (Human) - Antibody for Immunohistochemistry	100 µl
H-070-30	NEI-MCH (Prepro-MCH (131-165)) [Des-Asn7] (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-076-89	NERP-1 (Human) - Antibody for Immunohistochemistry	100 µl
H-076-91	NERP-2 (Human) - Antibody for Immunohistochemistry	100 µl
H-076-92	NERP-2 (Rat) - Antibody for Immunohistochemistry	100 µl
H-003-24	Nesfatin-1 (1-45) / Nesfatin-1 N-Terminal (Human) - Antibody for Immunohistochemistry	100 µl
H-003-22	Nesfatin-1 (1-82) (Rat) - Antibody for Immunohistochemistry	100 µl
H-003-25	Nesfatin-1 (46-82) (Human) - Antibody for Immunohistochemistry	100 µl
H-003-96	Nesfatin-1 C-Terminal (Human) - Antibody for Immunohistochemistry	100 µl
H-003-27	Nesfatin-3, C-Terminal (Human) - Antibody for Immunohistochemistry	100 µl
H-004-31	Neuroendocrine Protein 7B2, Pro (126-159) (Human, Porcine) - Antibody for Immunohistochemistry	50 µl
H-046-15	Neurokinin A / Neuromedin L / Substance K (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-046-26	Neurokinin B / Neuromedin K (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-045-88	Neuromedin S (Mouse) - Antibody for Immunohistochemistry	50 µl
H-045-87	Neuromedin S (Rat) - Antibody for Immunohistochemistry	100 µl
H-045-94	Neuromedin S, Prepro (70-103) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-045-93	Neuromedin S, Prepro (70-103) (Rat) - Antibody for Immunohistochemistry	100 µl
H-046-41	Neuromedin U (Rat) - Antibody for Immunohistochemistry	50 µl
H-045-96	Neuromedin U, prepro (104-136) (Human) - Antibody for Immunohistochemistry	100 µl
H-046-42	Neuromedin U-25 (Human) - Antibody for Immunohistochemistry	50 µl
H-046-40	Neuromedin U-25 (Porcine) - Antibody for Immunohistochemistry	50 µl
H-060-50	Neuronostatin-13 (Human, Porcine) - Antibody for Immunohistochemistry	100 µl
H-060-48	Neuronostatin-13 (Rat, Mouse) - Antibody for Immunohistochemistry	100 µl
H-060-53	Neuronostatin-19 (Human) - Antibody for Immunohistochemistry	100 µl
H-048-40	Neuropeptide AF (huNPAF) (Human) - Antibody for Immunohistochemistry	50 µl
H-005-51	Neuropeptide B-29 (NPB-29) (Human) - Antibody for Immunohistochemistry	50 µl
H-005-56	Neuropeptide B-29 (NPB-29) (Rat) - Antibody for Immunohistochemistry	50 µl
H-070-49	Neuropeptide EI (NEI) (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-005-81	Neuropeptide S (1-10) (Human) - Antibody for Immunohistochemistry	100 µl
H-005-77	Neuropeptide S, Prepro (23-67) (Human) - Antibody for Immunohistochemistry	100 µl
H-005-99	Neuropeptide S, Prepro (23-67) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-005-72	Neuropeptide S, Prepro (49-97) (Rat) - Antibody for Immunohistochemistry	100 µl

Catalog Number	Name	Size
H-005-61	Neuropeptide W-23 (Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-048-03	Neurotensin (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-048-02	Neurotensin, Prepro (24-55) (Human) Antibody for Immunohistochemistry	100 µl
H-021-58	Nociceptin (85-119), Pro / Orphanin FQ (85-119), Pro / Nocistatin-35 (Rat) - Antibody for Immunohistochemistry	50 µl
H-021-55	Nociceptin / Orphanin FQ (Human, Rat, Mouse, Ox) - Antibody for Immunohistochemistry	50 µl
NRS-IHC	Normal Rabbit Control Serum for Immunohistochemistry	500 µl
NRG-500	Normal Rabbit Purified IgG Control Antibody	500 µg
NRG-1000	Normal Rabbit Purified IgG Control Antibody	1 mg
H-002-62	NPDC-1 (137-176) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-64	NPDC-1 (54-81) (Human) - Antibody for Immunohistochemistry	100 µl
H-001-58	NPFF-1R (185-203) (Human) - Antibody for Immunohistochemistry	100 µl
H-001-60	NPFF-1R (414-430) (Human) - Antibody for Immunohistochemistry	50 µl
H-001-59	NPFF-2R (291-308) (Human) - Antibody for Immunohistochemistry	100 µl
H-049-03	NPY (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-049-07	NPY (Porcine) - Antibody for Immunohistochemistry	50 µl
H-031-92	Obestatin (Human, Monkey) - Antibody for Immunohistochemistry	100 µl
H-031-90	Obestatin (Rat, Mouse) - Antibody for Immunohistochemistry	100 µl
H-001-92	Olfactory Receptor 9A2 / GPR-139 (289-310) (Human) - Antibody for Immunohistochemistry	100 µl
H-003-36	Orexin A (16-33) Amide (Human, Rat, Mouse, Porcine, Ovine) - Antibody for Immunohistochemistry	50 µl
H-003-30	Orexin A / Hypocretin-1 (Human, Rat, Mouse, Porcine, Ovine, Bovine) - Antibody for Immunohistochemistry	50 µl
H-003-31	Orexin B / Hypocretin-2 (Human) - Antibody for Immunohistochemistry	50 µl
H-003-32	Orexin B / Hypocretin-2 (Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	50 µl
H-001-67	OT7T022 (243-264) (Neuropeptide FF Receptor 1) (Rat) - Antibody for Immunohistochemistry	50 µl
H-028-22	Oxyntomodulin (Human, Rat, Mouse) - Antibody for Immunohistochemistry	100 µl
H-028-16	Oxyntomodulin (Porcine) - Antibody for Immunohistochemistry	100 µl
H-051-01	Oxytocin (Human, Rat, Mouse, Bovine) - Antibody for Immunohistochemistry	50 µl
H-052-04	PACAP (6-27) (Human, Rat, Ovine) - Antibody for Immunohistochemistry	100 µl
H-052-02	PACAP 27 Amide (Human, Rat, Ovine) - Antibody for Immunohistochemistry	50 µl
H-052-08	PACAP 38 (28-38) (Human, Rat, Ovine) - Antibody for Immunohistochemistry	100 µl
H-052-05	PACAP 38 (Human, Rat, Ovine) - Antibody for Immunohistochemistry	50 µl
H-052-11	PACAP-Related Peptide (PRP) (Rat) - Antibody for Immunohistochemistry	50 µl
H-053-08	Pancreastatin (Porcine) - Antibody for Immunohistochemistry	50 µl
H-053-13	Pancreastatin / Chromogranin A (264-314)-Amide (Rat) - Antibody for Immunohistochemistry	50 µl
H-054-01	Pancreatic Polypeptide (Avian) - Antibody for Immunohistochemistry	50 µl
H-054-02	Pancreatic Polypeptide (Human) - Antibody for Immunohistochemistry	50 µl
H-054-07	Pancreatic Polypeptide (Rat) - Antibody for Immunohistochemistry	100 µl
H-001-88	PAR-3 (352-374) (Human) - Antibody for Immunohistochemistry	50 µl
H-055-08	Parathyroid Hormone (PTH) (1-34) (Human) - Antibody for Immunohistochemistry	50 µl
H-055-24	Parathyroid Hormone (PTH) (1-34) (Rat) - Antibody for Immunohistochemistry	50 µl
H-055-12	Parathyroid Hormone (PTH) (1-38) (Human) - Antibody for Immunohistochemistry	50 µl
H-055-13	Parathyroid Hormone (PTH) (1-44) (Human) - Antibody for Immunohistochemistry	50 µl
H-056-04	Parathyroid Hormone-related Protein (PTH-RP) (1-34) (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-056-01	Parathyroid Hormone-Related Protein (PTH-RP) (1-36) Amide [Tyr36] (Chicken) - Antibody for Immunohistochemistry	50 µl
H-024-58	Peptide F (Bovine) - Antibody for Immunohistochemistry	50 µl
H-064-01	PHI-27 (Porcine) - Antibody for Immunohistochemistry	50 µl
H-064-07	PHM-27 / PHI (Human) - Antibody for Immunohistochemistry	50 µl

Catalog Number	Name	Size
H-002-78	PK2L (48-68) (Human) - Antibody for Immunohistochemistry	100 µl
H-035-24	Preptin (Human) - Antibody for Immunohistochemistry	50 µl
H-074-01	Presenilin-1 (S182 Protein) (431-467) (Human) - Antibody for Immunohistochemistry	50 µl
H-002-35	Pro-Angiotensin-12 (Rat) - Antibody for Immunohistochemistry	100 µl
H-004-52	Pro-SAAS (221-242) / PEN (Human) - Antibody for Immunohistochemistry	100 µl
H-023-59	Prokineticin 1 (PK1) (81-105) [Cys(Acm)96] (Human) - Antibody for Immunohistochemistry	50 µl
H-008-50	Prolactin-Releasing Peptide-31 (PrRP-31) (Human) - Antibody for Immunohistochemistry	50 µl
H-008-52	Prolactin-Releasing Peptide-31 (PrRP-31) (Rat) - Antibody for Immunohistochemistry	50 µl
H-029-30	Proopiomelanocortin Precursor (POMC) (27-52) (Porcine) - Antibody for Immunohistochemistry	50 µl
H-023-69	Protein Inhibitor of nNOS (PIN) (24-58) (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-059-06	PYY (1-36)-Gly (Rat, Mouse, Porcine, Canine) - Antibody for Immunohistochemistry	100 µl
H-059-16	PYY (3-19) (Cys0) (Human) - Antibody for Immunohistochemistry	100 µl
H-059-17	PYY (3-19) (Cys20) (Human) - Antibody for Immunohistochemistry	100 µl
H-059-02	PYY (3-36) (Human) - Antibody for Immunohistochemistry	100 µl
H-059-04	PYY (3-36) (Rat, Mouse, Porcine, Canine) - Antibody for Immunohistochemistry	100 µl
H-059-09	PYY (3-36)-Gly (Canine, Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	100 µl
H-059-08	PYY (3-36)-Gly (Human) - Antibody for Immunohistochemistry	100 µl
H-059-01	PYY (Human) - Antibody for Immunohistochemistry	100 µl
H-059-03	PYY (Rat, Mouse, Porcine, Canine) - Antibody for Immunohistochemistry	100 µl
H-001-81	QRFP receptor AQ27 / SP9155 (401-433) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-002-74	RAR-related Orphan Receptor A isoform a (293-523) (Human) - Antibody for Immunohistochemistry	100 µl
H-035-62	Relaxin 2 (Human) - Antibody for Immunohistochemistry	100 µl
H-028-39	Resistin (43-65) (Human) - Antibody for Immunohistochemistry	50 µl
H-028-40	Resistin (51-108) (Human) - Antibody for Immunohistochemistry	50 µl
H-048-55	RFRP, Prepro (115-131) Amide (Human) - Antibody for Immunohistochemistry	50 µl
H-048-54	RFRP, Prepro (56-92) Amide (Human) - Antibody for Immunohistochemistry	100 µl
H-048-46	RFRP-3 (Human) - Antibody for Immunohistochemistry	100 µl
H-010-73	Salusin (50-mer) (Human) - Antibody for Immunohistochemistry	100 µl
H-010-67	Salusin-alpha (Human) - Antibody for Immunohistochemistry	100 µl
H-010-68	Salusin-beta (Human) - Antibody for Immunohistochemistry	100 µl
H-023-45	Sarafotoxin S6b (<i>Atractaspis engaddensis</i>) - Antibody for Immunohistochemistry	50 µl
H-023-46	Sarafotoxin S6c (<i>Atractaspis engaddensis</i>) - Antibody for Immunohistochemistry	50 µl
H-023-48	Sarafotoxin S6d (<i>Atractaspis engaddensis</i>) - Antibody for Immunohistochemistry	50 µl
H-006-10	Secreted Frizzled-related Protein 2 / SFRP 2 (25-50) (Human) - Antibody for Immunohistochemistry	100 µl
H-067-02	Secretin (Chicken) - Antibody for Immunohistochemistry	50 µl
H-067-03	Secretin (Human) - Antibody for Immunohistochemistry	50 µl
H-067-04	Secretin (Mouse) - Antibody for Immunohistochemistry	50 µl
H-067-05	Secretin (Porcine, Bovine) - Antibody for Immunohistochemistry	50 µl
H-067-06	Secretin (Rat) - Antibody for Immunohistochemistry	50 µl
H-047-95	Secretoneurin (Rat) - Antibody for Immunohistochemistry	50 µl
H-SEF-01	Sef (Human) - Antibody for Immunohistochemistry	100 µl
H-024-17	SNSR-1 (294-322) [Cys0] (Human) - Antibody for Immunohistochemistry	50 µl
H-060-03	Somatostatin (Human, Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	50 µl
H-060-14	Somatostatin-28 (Human, Rat, Mouse, Porcine) - Antibody for Immunohistochemistry	50 µl
H-006-22	SORL1 (2177-2214) (Human) - Antibody for Immunohistochemistry	100 µl
H-023-81	Spexin / NPQ (Human, Mouse, Bovine) - Antibody for Immunohistochemistry	100 µl
H-023-82	Spexin, Prepro (36-58) (Human) - Antibody for Immunohistochemistry	100 µl
H-023-83	Spexin, Prepro (73-116) (Human) - Antibody for Immunohistochemistry	100 µl
H-019-26	Stresscopin (SCP) (Human) - Antibody for Immunohistochemistry	100 µl

Catalog Number	Name	Size
H-019-27	Stresscopin Related Peptide (SRP) (Human) - Antibody for Immunohistochemistry	50 µl
H-019-27-P	Stresscopin Related Peptide (SRP) (Human) - Antibody for Immunohistochemistry (PPI Serum)	100 µl
H-061-05	Substance P (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-008-64	Synaptotagmin IV (73-97) (Human) - Antibody for Immunohistochemistry	100 µl
H-030-09	TGF-1, alpha (Rat) - Antibody for Immunohistochemistry	50 µl
H-002-06B	TGR-5 (298-330) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-06A	TGR-5 (298-330) (Human) - Antibody for Immunohistochemistry	50 µl
H-003-21	TRB3 (314-349) (Rat) - Antibody for Immunohistochemistry	100 µl
H-050-50	TRPM8 / CMR-1 (656-680) (Rat) - Antibody for Immunohistochemistry	50 µl
H-003-95	UCP3 (Uncoupling Protein 3) (110-124) (Human) - Antibody for Immunohistochemistry	50 µl
H-019-14	Urocortin (Human) - Antibody for Immunohistochemistry	50 µl
H-019-15	Urocortin (Rat) - Antibody for Immunohistochemistry	50 µl
H-019-30	Urocortin II (Human) - Antibody for Immunohistochemistry	50 µl
H-019-24	Urocortin II (Mouse) - Antibody for Immunohistochemistry	50 µl
H-019-28	Urocortin III (Human) - Antibody for Immunohistochemistry	50 µl
H-019-29	Urocortin III (Mouse) - Antibody for Immunohistochemistry	100 µl
H-005-18	Urodilatin / ANP (95-126) - Antibody for Immunohistochemistry	50 µl
H-071-05	Urotensin II (Human) - Antibody for Immunohistochemistry	50 µl
H-071-08	Urotensin II (Mouse) - Antibody for Immunohistochemistry	50 µl
H-071-09	Urotensin II (Rat) - Antibody for Immunohistochemistry	50 µl
H-071-17	Urotensin II Related Peptide (Human, Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-071-24	Urotensin II, Prepro (21-40) (Human) - Antibody for Immunohistochemistry	100 µl
H-071-23	Urotensin II, Prepro (41-68) (Human) - Antibody for Immunohistochemistry	100 µl
H-071-22	Urotensin II, Prepro (71-84) (Human) - Antibody for Immunohistochemistry	100 µl
H-V007	VACHt (Human, Rat, Mouse, Monkey) - Antibody for Immunohistochemistry	50 µl
H-V006	VACHt (Rat, Mouse) - Antibody for Immunohistochemistry	50 µl
H-V005	VACHt C-terminal Peptide (Human) - Antibody for Immunohistochemistry	50 µl
H-002-95	VASH 1 Protein (177-204) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-96	Vasohibin [Cys0] (286-299) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-94	Vasohibin [Cys0] (336-365) (Human) - Antibody for Immunohistochemistry	100 µl
H-065-07	Vasopressin (AVP) [Arg8] (Human, Rat, Mouse, Ovine) - Antibody for Immunohistochemistry	50 µl
H-053-34	Vasostatin I (17-43) (Human) - Antibody for Immunohistochemistry	100 µl
H-053-35	Vasostatin I (17-76) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-66	Vaspin (386-414) (Human) - Antibody for Immunohistochemistry	100 µl
H-002-67	Vaspin (63-104) (Human) - Antibody for Immunohistochemistry	100 µl
H-064-16	VIP (Human, Rat, Mouse, Porcine, Ovine) - Antibody for Immunohistochemistry	50 µl
H-003-93	Visfatin (397-422) (Human) - Antibody for Immunohistochemistry	100 µl
H-003-84B	Visfatin (400-450) (Human) - Antibody for Immunohistochemistry	200 µl
H-003-84	Visfatin (400-450) (Human) - Antibody for Immunohistochemistry	100 µl
H-003-82	Visfatin (412-431) (Human) - Antibody for Immunohistochemistry	100 µl
H-003-79	Visfatin (467-491) (Human) - Antibody for Immunohistochemistry	100 µl
H-003-91	Visfatin / PBEF Isoform I, His (Rec.) (Mouse) - Antibody for Immunohistochemistry	100 µl
H-V001	VMAT 1 (Human) - Antibody for Immunohistochemistry	50 µl
H-V002	VMAT 1 (Rat) - Antibody for Immunohistochemistry	50 µl
H-V003	VMAT 2 (Human) - Antibody for Immunohistochemistry	50 µl
H-V008	VMAT 2 (Mouse) - Antibody for Immunohistochemistry	50 µl
H-V004	VMAT 2 (Rat) - Antibody for Immunohistochemistry	50 µl
H-046-74	Xenin 25 (Human, Rat, Mouse) - Antibody for Immunohistochemistry	100 µl