

Phoenix Custom Peptides

SYNTHESIZING THE IMPOSSIBLE

At Phoenix, we have over 40 years of experience synthesizing diverse and difficult peptides from insoluble tripeptides to those with multiple disulfide bridges and linked a- and b- chains

ABOUT OUR SERVICES

Phoenix's innovative, widely-cited peptide ligands are always the first to be produced commercially for the research market.


Phosphorylation, sulfonation, octanoylation, unusual amino acids, etc.

HPLC and Mass Spec standard


A few of our achievements:

Ghrelin: A 28-mer peptide with n-Octanoyl-Ser at the third position.

AgRP (83-132)-NH₂: A 50-mer polypeptide with 5 disulfide bridges



AGRP (83-132) NH₂



GHRELIN (HUMAN)

Send quote request to:
quotes@phoenixpeptide.com



PHOENIX PHARMACEUTICALS, INC.
330 BEACH ROAD, BURLINGAME CA, 94010, USA
PHONE: (650) 558-8898 EMAIL: info@phoenixpeptide.com
WWW.PHOENIXPEPTIDE.COM

PHOENIX EUROPE GMBH
VIKTORIASTRASSE 3-5, D-76133 KARLSRUHE, GERMANY
PHONE: +49-721-1611950 EMAIL: germany@phoenixpeptide.com
WWW.PHOENIXPEPTIDE.COM