

# The Peak of Pipetting

## New! 300 uL XL Multi-Volume Pipette Tip

- 2-300 uL—non-filtered volume range
- 2-200 uL—filtered volume range

#### One tip does it all!

- Graduations at 10, 50, 100, 150 and 200  $\mu$ L—eliminate the need to keep multiple tips in your inventory.
- Universal fit for broad compatibility with most major pipettor brands.
- Filter placement is optimized to accommodate 20  $\mu$ L, 100  $\mu$ L, and 200  $\mu$ L pipettors; non-filtered tips are also compatible with 300  $\mu$ L pipettors.

### **Complete sample recovery!**

- 91 mm (3.58 in) in length—easily reach the bottom of 12 x 75 mm tubes without risk of the pipettor nose cone touching the interior of the tube.
- Manufactured from virgin polypropylene to ensure purity and provide excellent optical clarity.
- Vertex NoStick® resin technology creates a low-binding surface for maximum sample recovery.
- Certified free of DNAse, RNAse, DNA and PCR inhibitors.
- Sterile product is also pyrogen-free.
- Available in standard 8 x 12 SBS-format rack.
- Patented press-to-open hinged lid enables easy, one-hand operation.

#### **Vertex Multi-Volume Pipette Tips**

Item #	Max Volume	Description	Color	Packaging
4280N00S	300 μL	300 μL XL Bulk	Natural	1000/pack 10 packs/case
4287N00S	300 μL	300 μL XL Racked	Natural	2 racks of 96/unit 5 units/case
4287NS0S	300 µL	300 μL XL Racked, Sterile	Natural	2 racks of 96/unit 5 units/case
4287NSFS	200 μL	200 μL XL Filtered Racked, Sterile	Natural	2 racks of 96/unit 5 units/case

