

Open Biosystems 社 Clone Collection / RNAi Collection 価格表

(2011年3月11日現在)

メーカーコード	販売用品名	新定価
BHS1023	GenMap Human BAC Clone	14,000
BHS1033	CCAP Human BAC Clone	14,000
BMM1035	Mere Mouse BAC Clone	14,000
EBT3996	Bos taurus IMAGE Clone	12,000
EDM1133	Drosophila Gene Collection Clone	17,000
EDR1052	IMAGE Zebrafish cDNA Clone	12,000
EDR4009	NEI Zebrafish Clone	17,000
EDR4422	Zebrafish EST cDNA (Exelixis)	17,000
EDR5649	Zebrafish ZGC PCR-Rescued cDNA	17,000
EFR1003	IMAGE Fugu cDNA Clone	12,000
EHS1001	IMAGE Human cDNA Clone	12,000
EHS1019	Human Cystic Fibrosis Clone	14,000
EHS1031	Human Eye cDNA Clone	14,000
EHS1193	NEI Bank Clone	17,000
EHS1194	UI Human Placenta EST Clone	17,000
EMM1002	IMAGE Mouse cDNA Clone	12,000
EMM1005	IMAGE Mouse Lymphochip Clone	12,000
EMM1032	BMAP Non-Redundant Clone	14,000
EMM1212	Organ of Corti Clone	17,000
EMM1384	Mouse Eye cDNA Clone	17,000
EMM3836	NIA Mouse Clone	17,000
EOC4011	cDNA: NEI rabbit clone	17,000
ERM1145	IMAGE Rhesus Monkey cDNA Clone	12,000
ERN1017	Rat Non-Redundant Clone	17,000
ERN1146	IMAGE Rat cDNA Clone	12,000
ERN1184	IMAGE Rat EST Clone	12,000
EUI1044	cDNA EST Clone	17,000
EXL1051	IMAGE Xenopus cDNA Clone	12,000
EXT1168	Xenopus Tropicalis IMAGE Clone	12,000
ICN1378	Incyte Dog cDNA Clone	18,000
IHS1380	Incyte FL Human cDNA Clone	18,000
IHS1382	Incyte Easy-to-Spot Human cDNA Clone	18,000
IRM1377	Incyte Macaque cDNA Clone	18,000
IRN1383	Incyte Rat cDNA Clone	18,000
MBT4353	Cow MGC Verified FL cDNA	16,000
MDR1734	ZGC Clone	17,000
MDR1738	IMCB Zebrafish Clone	16,000
MHS1010	MGC Human Verified FL Clone	17,000
MHS1011	MGC Human Verified FL Clone	17,000
MHS1768	Human MGC PCR-Rescued cDNA	17,000
MHS4426	Human MGC Verified FL cDNA clone (IRCM)	17,000
MHS4768	Human MGC verified FL cDNA (IRQL)	17,000
MHS4771	Human MGC verified FL cDNA (IRQM)	17,000
MMM1013	MGC Mouse Verified FL Clone	17,000
MMM4769	Mouse MGC verified FL cDNA (IRQL)	17,000
MRN1768	MGC Rat Verified FL Clone	17,000
MRN4770	Rat MGC verified FL cDNA (IRQL)	17,000
MXL1736	Xenopus laevis Verified FL Clone	17,000
MXT1765	Xenopus tropicalis Verified FL Clone	17,000
OCE1182	C. elegans ORF Clone	14,000
OCE4818	C. elegans Transcription Factor Clone - E. coli	14,000
ODR5433	Zebrafish ORFeome Collaboration Clone	69,000
OEC4626	E. coli SPA-Tagged Strain	14,000
OEC4629	E. coli TAP-Tagged Strain	14,000
OHS1770	Human ORF Clone	52,000
OHS4464	Human Synthesized Full-Length ORF(MGC - IRCH)	17,000

メーカーコード	販売用品名	新定価
OHS4766	Human synthesized ORF clone (IRCV)	17,000
OHS4903	Human SH2 Domain Clone	14,000
OHS5834	Precision LentiORF Expression Transfection Starter Kit	207,000
OHS5835	Precision LentiORF Expression Transduction Starter Kit	311,000
OHS5836	Precision LentiORF Expression Viral Particle Starter Kit	233,000
OHS5893	Human ORFeome Collaboration Clone (with Native Stop Codon),Glycerol	69,000
OHS5894	Human ORFeome Collaboration Clone (w/o Stop Codon),Glycerol	69,000
OHS5897	Precision LentiORF Individual Clone (with Native Stop Codon),Glycerol	78,000
OHS5898	Precision LentiORF Individual Clone (w/o Stop Codon),Glycerol	78,000
OHS5899	Precision LentiORF Individual Clone (with Native Stop Codon),Viral Particles	130,000
OHS5900	Precision LentiORF individual clone (w/o Stop Codon),Viral Prticles	130,000
OMM4767	Mouse synthesized ORF clone (IRCW)	17,000
OMM5895	Mouse ORFeome Collaboration Clone (with Native Stop Codon),Glycerol	69,000
OMM5896	Mouse ORFeome Collaboration Clone (w/o Stop Codon),Glycerol	69,000
PCE1182	C.elegans Promoters v1.1 Clone	14,000
PEC3876	E. coli Promoter Clone	14,000
RCE1182	C.elegans RNAi Clone	19,000
RHS1764	Human retroviral shRNAmir individual clone	37,000
RHS3600	pSM2 Retroviral shRNAmir Starter Kit-Luciferase,RNAintro	156,000
RHS3601	pSM2 Retroviral shRNAmir Starter Kit-eGFP,RNAintro	156,000
RHS3979	TRC Human shRNA Construct	36,000
RHS4287	GIPZ Lentiviral shRNAmir Starter Kit,RNAintro	156,000
RHS4430	Human pGIPZ shRNAmir individual clone	37,000
RHS4529	Human pSM2 retroviral shRNAmir target gene set	75,000
RHS4531	Human GIPZ lentiviral shRNAmir target gene set	75,000
RHS4533	Human pLKO.1 lentiviral shRNAmir targetgene set	56,000
RHS4696	Human TRIPZ inducible lentiviral shRNAmir individual clone	37,000
RHS4740	Human TRIPZ lentiviral inducible shRNAmir target gene set	75,000
RHS4741	TRIPZ Lentiviral Inducible shRNAmir Starter Kit,RNAintro	156,000
RHS4917	Human pSMP retroviral shRNAmir individual clone	37,000
RHS5043	pSMP Retroviral shRNAmir Starter Kit-Luciferase,RNAintro	156,000
RHS5086	GIPZ Lentiviral shRNAmir Transduction Kit,RNAintro,for Academic	232,000
RHS5087	TRIPZ Lentiviral Inducible shRNAmir Transduction Kit,RNAintro	232,000
RMM1766	Mouse retroviral shRNAmir individual clone	37,000
RMM3981	TRC Mouse shRNA Construct	36,000
RMM4431	Mouse pGIPZ shRNAmir individual clone	37,000
RMM4530	Mouse pSM2 retroviral target gene shRNAmir set	75,000
RMM4532	Mouse GIPZ lentiviral target gene shRNAmir set	75,000
RMM4534	Mouse pLKO.1 lentiviral target gene shRNA set	56,000
VGH5518	Lentiviral shRNAmir,Human,Viral Particle,Individual Clone,GIPZ	135,000
VGM5520	Lentiviral shRNAmir,Mouse,Viral Particle,Individual Clone,GIPZ	135,000
YSC1021	Yeast Knock Out Clone	10,000
YSC1070	HA Tag Clone	10,000
YSC1073	Mini- Transposon Clone	10,000
YSC1178	TAP Tagged ORF Clone	10,000
YSC1180	Yeast Tet-promoters Hughes (yTHC) Strain	10,000
YSC1202	Yeast mTn Individual Plasmid	10,000
YSC1210	Yeast Tet-promoters Hughes Parental Strain	10,000
YSC1764	HA Tag Parental Strain	10,000
YSC3867	Yeast ORF Clone - E.coli	10,000
YSC3869	Yeast ORF Clone	10,000
YSC4034	Yeast Magic Marker Strain (KanMX)	10,000
YSC4427	Yeast Magic Marker Strain (URA3)	10,000
YSC4515	Yeast GST-Tagged Strain	10,000
YSC5093	Yeast DAmP individual strain - diploid	10,000
YSC5094	Yeast DAmP individual strain - haploid	10,000
YSC5095	Yeast Cross & Capture strain - Bait MATa	10,000
YSC5096	Yeast Cross & Capture strain - Prey MATalpha	10,000
YSC5097	Yeast Protein Interactome - Individual Bait F[1,2]	10,000
YSC5098	Yeast Protein Interactome - Individual Prey F[3]	10,000