

Anti-Zebrafish HSP70 : ZHSP67-42 (Mouse monoclonal antibody)

Code : ZHSP6742
Lot :
Lot size : 1 mL
Isotype : IgG
Antigen : HSP 70 (Zebrafish)
Purity : >80% by SDS-PAGE
Concentration : 0.1 mg /mL
Appearance : Clear liquid
Storage : Short term, 4°C
Long term, -20°C
Avoid freeze thaw cycles
Preservative : 0.1% Sodium azide
Presentation : Supplied as a liquid in 0.01M Tris,0.15M NaCl,
0.1% Sodium azide, pH:7.5

EnBioTec Laboratories, Co., Ltd.

Chiyoda-parion Bldg.,6th Floor, 2-3-16 Kandasuda-chou, Chiyoda-ku, Tokyo 101-0041, Japan

Tel:03-5297-7155 Fax:03-5297-0581 E-mail: support@enbiotec.co.jp