

nippi Low Endotoxin Gelatins

★★★Line Up★★★

Products	Type	Feature	Viscosity	Gel Formation
MediGelatin	HMG-BP	High Gel Strength (over 250bloom)	High	Gel
Nippi High Grade Gelatin	APAT	Middle Molecular Weight (about 60,000)	↕	↕
	AP	Low Molecular Weight (about 8,000)		




- ★ From Porcine Skin
- ★ Sterilized
- ★ Low Endotoxin
(Less Than 1 O E U / g)


Gelatin which is extremely reduced endotoxin content and high gel strength for medical use.

Therefore “nippi Low Endotoxin Gelatins” inactivated for immune system macrophage that reacted to endotoxin.

Used for : Stabilizer for injection , Drug delivery systems ,
Tissue engineering etc...



NO release by murine macrophages activated with LPS or gelatins



Contact us

Address :

**1-1-1 Senju-midoricho , Adachi-ku
Tokyo 120-8601 , Japan
Nippi the Protein Engineering Office**

Phone & Fax : +81(3)3888-5180

E-mail : pex@nippi-inc.co.jp