

## ARG1 Antibody (C-term)

Peptide Affinity Purified Rabbit Polyclonal Antibody (Pab)

Catalog # AP8976b

### Specification

#### ARG1 Antibody (C-term) - Product Information

Application	WB, IHC-P, FC,E
Primary Accession	<a href="#">P05089</a>
Reactivity	Human
Host	Rabbit
Clonality	Polyclonal
Isotype	Rabbit Ig
Clone Names	RB23044
Antigen Region	293-322

#### ARG1 Antibody (C-term) - Additional Information

Gene ID 383

#### Other Names

Arginase-1, Liver-type arginase, Type I arginase, ARG1

#### Target/Specificity

This ARG1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 293-322 amino acids from the C-terminal region of human ARG1.

#### Dilution

WB~~1:1000  
IHC-P~~1:50~100  
FC~~1:10~50

#### Format

Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

#### Storage

Maintain refrigerated at 2-8°C for up to 6 months. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.


#### Precautions

ARG1 Antibody (C-term) is for research use only and not for use in diagnostic or therapeutic procedures.


#### ARG1 Antibody (C-term) - Protein Information

Name ARG1

Cellular Location



Western blot analysis of ARG1 Antibody (C-term) (Cat. #AP8976b) in MDA-MB231 cell line lysates (35ug/lane). ARG1 (arrow) was detected using the purified Pab.




Formalin-fixed and paraffin-embedded human hepatocarcinoma reacted with ARG1 Antibody (C-term), which was peroxidase-conjugated to the secondary antibody, followed by DAB staining. This data demonstrates the use of this antibody for immunohistochemistry; clinical relevance has not been evaluated.

Cytoplasm.

### ARG1 Antibody (C-term) - Protocols

Provided below are standard protocols that you may find useful for product applications.

- [Western Blot](#)
- [Blocking Peptides](#)
- [Dot Blot](#)
- [Immunohistochemistry](#)
- [Immunofluorescence](#)
- [Immunoprecipitation](#)
- [Flow Cytometry](#)
- [Cell Culture](#)



ARG1 Antibody (C-term) (Cat. #AP8976b) flow cytometric analysis of MDA-MB231 cells (right histogram) compared to a negative control cell (left histogram). FITC-conjugated goat-anti-rabbit secondary antibodies were used for the analysis.

### ARG1 Antibody (C-term) - Background

Arginase catalyzes the hydrolysis of arginine to ornithine and urea. At least two isoforms of mammalian arginase exist (types I and II) which differ in their tissue distribution, subcellular localization, immunologic crossreactivity and physiologic function. The type I isoform encoded by this gene, is a cytosolic enzyme and expressed predominantly in the liver as a component of the urea cycle. Inherited deficiency of this enzyme results in argininemia, an autosomal recessive disorder characterized by hyperammonemia.

### ARG1 Antibody (C-term) - References

Haraguchi, Y., et al., J. Clin. Invest. 86 (1), 347-350 (1990)