

Biotin Single-Shot Antibody Labeling Kit

Catalog #	Source:	Lot #	Size:
KT1001a			
Concentration:	Accession:	Clone Name:	Isotype:
Applications:		Reactivity:	MW (kDa):

Target/Specificity:

Other Names:

Application Data:



This product is offered at a low price as part of the Abgent Product Investigation Program. Additionally, submit your positive data for this or any Abgent antibody to our [Data Rewards Program](#) and receive a free gift.

Background:

Abgent's Biotin Single-Shot Antibody Labeling Kit offers the complete solution to the quantitative biotinylation of 100 ug of antibody . The kit incorporates traceable biotin and allows for direct spectrophotometric verification of biotin incorporation. A simple λ 280 and λ 354 measurement determines the number of biotins per antibody molecule.

[User Manual](#)

Background References:

Product Citations:

Application Notes:

Format:

Storage:

Precautions:

This product is for research use only. Not for use in diagnostic or therapeutic procedures.