

Datasheet

RCC1 monoclonal antibody (M01), clone 2F1

Catalog Number: H00001104-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant RCC1.

Clone Name: 2F1

Immunogen: RCC1 (AAH07300, 312 a.a. ~ 421 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

EGKAYSLGRAEYGRGLGEGAEKSIPTLISRLPAVSS
VACGASVGYAVTKDGRVFAWGMGTNYQLGTGQDED
AWSPVEMMGKQLENRVVLSVSSGGQHTVLLVKDKEQ
S

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, S-ELISA, WB-Ce, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 1104

Gene Symbol: RCC1

Gene Alias: CHC1, RCC1-I