

Datasheet

FOXC2 monoclonal antibody (M02A), clone 2H3

however, it may play a role in the development of mesenchymal tissues. [provided by RefSeq]

Catalog Number: H00002303-M02A

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant FOXC2.

Clone Name: 2H3

Immunogen: FOXC2 (NP_005242.1, 421 a.a. ~ 501 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

AASWYLNHSGDLNHLPGHTFAAQQTFFPNVREMFNS
HRLGIENSTLGESQVSGNASCQLPYRSTPPLYRHAAP
YSYDCTKY

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgM Lambda

Storage Buffer: In ascites fluid

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2303

Gene Symbol: FOXC2

Gene Alias: FKHL14, LD, MFH-1, MFH1

Gene Summary: This gene belongs to the forkhead family of transcription factors which is characterized by a distinct DNA-binding forkhead domain. The specific function of this gene has not yet been determined;