

Datasheet

GPR3 monoclonal antibody (M01), clone 3B4-G3

Gene Alias: ACCA

Catalog Number: H00002827-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant GPR3.

Clone Name: 3B4-G3

Immunogen: GPR3 (AAH32702, 1 a.a. ~ 330 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

```
MMWGAGSPLAWLSAGSGNVNVSSVGP AEGPTGPAA  
PLPSPKAWDVVLCISGTLVSCENALVVAIIVGTPAFRAP  
MFLLVGSLAVADLLAGLGLVLHFAAVFCIGSAEMSLVL  
VGLVMAFTASIGSLLAITVDRYLSLYNALTYSETTVT  
RTYVMLALVWGGALGLLLPVLAWNCLDGLTTCGVVY  
PLSKNHLVVLAIAFFMVF GIMLQLYAQICRIVCRHAQQI  
ALQRHLLPASHYVATRKG IATLAVVLGAF AACWLPFTV  
YCLLGDAHSPPLYTYL TLLPATYNSMINPIIYAFRNQDV  
QKVLWAVCCCCSSSKIPFRSRSPSDV
```

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, RNAi-Ab, S-ELISA, WB-Ce, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 2827

Gene Symbol: GPR3