

Datasheet

HLA-C purified MaxPab mouse polyclonal antibody (B01P)

Catalog Number: H00003107-B01P

Regulatory Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human HLA-C protein.

Immunogen: HLA-C (AAH02463, 1 a.a. ~ 366 a.a) full-length human protein.

Sequence:

```
MRVMAPRTLILLLSGALALTETWACSHSMRYFYTAVS  
RPGRGEPFRFIAVGYVDDTQFVRFSDAASPRGEPRAP  
WVEQEGPEYWDRETQKYKRQAQTDRVSLRNLRGYY  
NQSEAGSHTLQWMYGCDLGPDRLLRGYDQSAYDG  
KDYIALNEDLRSWTAADTAAQITQRKWEAARAAEQQR  
AYLEGTCEVWLRRLRYLENGKETLQRAEHPKTHVTHHLV  
SDHEATLRCWALGFYPAEITLTWQRDGEDQTQDTELV  
ETRPAGDGTQKWAAVVVPSGEEQRYTCHVQHEGLP  
EPLTLRWEPSQPTIPIVGIVAGLAVLAVLAVLGAVVAV  
VMCRRKSSGGKGGSCSQAASSNSAQQSDESLIACKA
```

Host: Mouse

Reactivity: Human

Applications: WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 3107

Gene Symbol: HLA-C

Gene Alias: D6S204, FLJ27082, HLA-Cw, HLA-Cw12, HLA-JY3, HLC-C, PSORS1

Gene Summary: HLA-C belongs to the HLA class I

heavy chain paralogues. This class I molecule is a heterodimer consisting of a heavy chain and a light chain (beta-2 microglobulin). The heavy chain is anchored in the membrane. Class I molecules play a central role in the immune system by presenting peptides derived from endoplasmic reticulum lumen. They are expressed in nearly all cells. The heavy chain is approximately 45 kDa and its gene contains 8 exons. Exon one encodes the leader peptide, exons 2 and 3 encode the alpha1 and alpha2 domain, which both bind the peptide, exon 4 encodes the alpha3 domain, exon 5 encodes the transmembrane region, and exons 6 and 7 encode the cytoplasmic tail. Polymorphisms within exon 2 and exon 3 are responsible for the peptide binding specificity of each class one molecule. Typing for these polymorphisms is routinely done for bone marrow and kidney transplantation. Over one hundred HLA-C alleles have been described [provided by RefSeq]