

Datasheet

INSRR monoclonal antibody (M01), clone 6E6

Catalog Number: H00003645-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant INSRR.

Clone Name: 6E6

Immunogen: INSRR (NP_055030, 651 a.a. ~ 760 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

LYLNDYCHRGLRLPTSNNDRFDGEDGDPEAEMESD
CCPCQHPPPGQVLPPLEAQEASFQKKFENFLHNAITIP
ISPWKVTINKSPQRDSGRHRRRAAGPLRLGGNSSDF

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 3645

Gene Symbol: INSRR

Gene Alias: IRR