

Datasheet

PGM3 monoclonal antibody (M01A), clone S3

Catalog Number: H00005238-M01A

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant PGM3.

Clone Name: S3

Immunogen: PGM3 (AAH01258.1, 1 a.a. ~ 542 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MDLGAITKYSALHAKPNGLILQYGTAGFRTKAEHLDHV
MFRMGLLAVLRSKQTKSTIGVMVTASHNPEEDNGVKL
VDPLGEMLAPSWEEHATCLANAEEQDMQRVLIDISEK
EAVNLQQDAFVVIGRDTRPSSEKLSQSVIDGVTVLGG
QFHDYGLLTTPQLHYMVYCRNTGGRYGKATIEGYYQK
LSKAFVELTKQASCSGDEYRSLKVDCANGIGALKLRE
MEHYFSQGLSVQLFNDGSKGKLNHLGADFKSHQK
PPQGMEIKSNERCCSFDGDADRIVVYHDADGHFHLI
DGDKIATLISSFLKELLVEIGESLNIGVVQTAYANGSSTR
YLEEVMKVPVYCTKTGVKHLHHAQEFDIGVYFEANG
HGTALFSTAVEMKIKQSAEQLEDKCRKAAMLENIIDLF
NQAAGDAISDMLVIEAILALKGLTVQQWDALYTDLPNR
QLKVQVADRRVISTTNAERQAVTPPGLQEAINDLVKKY
KLSRAFVRPSGTEDVVRVYAEADSQESADHLAHEVSL
AVFQLAGGIGERPQPGF

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In ascites fluid

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5238

Gene Symbol: PGM3

Gene Alias: AGM1, DKFZp434B187, PAGM