

Datasheet

PTN monoclonal antibody (M02), clone 2E3

Catalog Number: H00005764-M02

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant PTN.

Clone Name: 2E3

Immunogen: PTN (NP_002816, 45 a.a. ~ 154 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

SDCGEWQWSVCVPTSGDCGLGTREGTRTGAECKQT
MKTQRCKIPCNWKKQFGAECKYQFQAWGECDLNTAL
KTRTGSLKRALHNAECQKTVTISKPCGKLTKPKQAES
K

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 5764

Gene Symbol: PTN

Gene Alias: HARP, HBGF8, HBNF, NEGF1