

Datasheet

RPA3 monoclonal antibody (M01), clone 1F4

Catalog Number: H00006119-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant RPA3.

Clone Name: 1F4

Immunogen: RPA3 (AAH05264, 12 a.a. ~ 121 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

INAGMLAQFIDKPVCFVGRLEKIHPTGKMFILSDGEGK
NGTIELMEPLDEEISGIVEVVGRVTAKATILCTSYYVQFK
EDSHPFDLGLYNEAVKIIHDFPQFYPLGIVQHD

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, RNAi-Ab, S-ELISA, WB-Ce, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 lambda

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 6119

Gene Symbol: RPA3

Gene Alias: REPA3