

Datasheet

SH3BGR monoclonal antibody (M01), clone 3B7

Catalog Number: H00006450-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant SH3BGR.

Clone Name: 3B7

Immunogen: SH3BGR (AAH06371, 1 a.a. ~ 239 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

```
MPLLLLGETEPLKLERDCRSPVEPWAAASPDALACL  
CHCQDLSSGAFPNRGVLLGGVLFPTVEMVIKVFVATSS  
GSAIRKKQQE VVGFLEANKIDFKELDIAGDEDNRRWM  
RENVPGKPKQNGIPLPPQIFNEEQYCGDFDSFFSAK  
EENIIYSFLGLAPPPDSKGSEKAEEGGETEAQKEGSED  
VGNLPEAQEKNEEEGETATEETEEIAMEGAEGEAE  
EETAEGEEPGEDEDS
```

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 6450

Gene Symbol: SH3BGR

Gene Alias: 21-GARP