

## Datasheet

### SHC1 monoclonal antibody (M01), clone 3F4

**Catalog Number:** H00006464-M01

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a partial recombinant SHC1.

**Clone Name:** 3F4

**Immunogen:** SHC1 (AAH33925, 171 a.a. ~ 280 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

QTPSHLGATLPVGGQVGGDPEVRKQMPPPPPCGRE  
LFDDPSYVNVQNLDKARQAVGGAGPPNPAINGSAPR  
DLFDMKPFEDALRVPPPPQSVMAEQLRGEPWFHGK  
LS

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, IF, IHC-P, PLA-Ce, S-ELISA, WB-Ce, WB-Re  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG1 Kappa

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 6464

**Gene Symbol:** SHC1

**Gene Alias:** FLJ26504, SHC, SHCA