

Datasheet

SNRPA monoclonal antibody (M01), clone 3F9-1F7

Gene Alias: U1A

Catalog Number: H00006626-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant SNRPA.

Clone Name: 3F9-1F7

Immunogen: SNRPA (AAH00405, 1 a.a. ~ 282 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MAVPETRPNHTIYINNLLNEKIKKDELKKSLEYAIFSQFGQI
LDILVSRSLKMRGQAFVIFKEVSSATNALRSMQGFPPY
DKPMRIQYAKTDSIIAKMKGTFVERDRKREKRKPKSQ
ETPATKKAVQGGGATPVVGAVQGPVPGMPPMTQAPR
IMHHMPGQPPYMPPPGMIPPPGLAPGQIPPGAMPPQ
QLMPGQMPPAQPLSENPPNHILFTNLPEETNELMLS
MLFNQFPGFKEVRLVPRHDIAFVEFDNEVQAGAARD
ALQGFKITQNNAMKISFAKK

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, S-ELISA, WB-Ce,
WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 6626

Gene Symbol: SNRPA