

## Datasheet

### SP1 monoclonal antibody (M01A), clone 4H6

**Catalog Number:** H00006667-M01A

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a partial recombinant SP1.

**Clone Name:** 4H6

**Immunogen:** SP1 (NP\_612482, 89 a.a. ~ 198 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

GTGELDLTATQLSQGANGWQIISSSSGATPTSKEQSG  
SSTNGSNGSESSKNRTVSGGQYVVAAPNLQNNQVL  
TGLPGVMPNIQYQVIPQFQTVGQQLQFAATGAQVQ  
Q

**Host:** Mouse

**Reactivity:** Human, Mouse

**Applications:** ELISA, WB-Ce, WB-Re  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG

**Storage Buffer:** In ascites fluid

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 6667

**Gene Symbol:** SP1

**Gene Alias:** -