

Datasheet

TFRC monoclonal antibody (M01A), clone 1E6

Catalog Number: H00007037-M01A

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant TFRC.

Clone Name: 1E6

Immunogen: TFRC (AAH01188, 68 a.a. ~ 168 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

YGTIAVIVFFLIGFMIGYLG YCKGVEPKTECERLAGTES
PVREEPGEDFPAARRLYWDDLKRLSEKLDSTDFGT
IKLLNENSYPREAGSQDENLALY

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2b Kappa

Storage Buffer: In ascites fluid

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 7037

Gene Symbol: TFRC

Gene Alias: CD71, TFR, TFR1, TRFR