

Datasheet

THY1 polyclonal antibody (A01)

Catalog Number: H00007070-A01

Regulatory Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a partial recombinant THY1.

Immunogen: THY1 (NP_006279, 26 a.a. ~ 134 a.a) partial recombinant protein with GST tag.

Sequence:

TACLVDQSLRLDCRHENTSSSPIQYEFSLTRETKKHVL
FGTVGVPEHTYRSRTNFTSKYNMKVLYLSAFTSKDEG
TYTCALHHSRGHSPPISSQNVTVLRDKLVKCEGIS

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: 50 % glycerol

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 7070

Gene Symbol: THY1

Gene Alias: CD90, FLJ33325