

## Datasheet

### THY1 purified MaxPab mouse polyclonal antibody (B02P)

**Catalog Number:** H00007070-B02P

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse polyclonal antibody raised against a full-length human THY1 protein.

**Immunogen:** THY1 (ABM85939.1, 1 a.a. ~ 161 a.a) full-length human protein.

**Sequence:**

MNLAISIALLLTVLQVSRGQKVTSLTACLVDQSLRLDCR  
HENTSSSPIQYEFSLTRETKKHVLFGTVGVPEHTYRSR  
TNFTSKYNMKVLVLSAFTSKDEGTYTCALHHSGHSPPI  
SSQNVTVLRDKLVKCEGISLLAQNTSWLLLLLSLSLLQ  
ATDFMSL

**Host:** Mouse

**Reactivity:** Human

**Applications:** WB-Tr

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 7070

**Gene Symbol:** THY1

**Gene Alias:** CD90, FLJ33325