

Datasheet

SF1 monoclonal antibody (M01A), clone 2E12

A. Nucleic Acids Res. 2015 Sep 29.

Catalog Number: H00007536-M01A**Regulatory Status:** For research use only (RUO)**Product Description:** Mouse monoclonal antibody raised against a partial recombinant SF1.**Clone Name:** 2E12**Immunogen:** SF1 (NP_004621, 1 a.a. ~ 110 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.**Sequence:**MATGANATPLDFPSKRRKRSRWNQDTMEQKTVIPGM
PTVIPGLTREQERAYIVQLQIEDLTRKLRTGDLGIPPN
PEDRSPSPEPIYNSEGKRLNTRERTRKKLEERH**Host:** Mouse**Reactivity:** Human, Mouse**Applications:** ELISA, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols**Isotype:** IgG2a Kappa**Storage Buffer:** In ascites fluid**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.**Entrez GeneID:** 7536**Gene Symbol:** SF1**Gene Alias:** D11S636, ZFM1, ZNF162**References:**

1. Mammalian splicing factor SF1 interacts with SURP domains of U2 snRNP-associated proteins. Crisci A, Raleff F, Bagdiul I, Raabe M, Urlaub H, Rain JC, Kramer