

Datasheet

BRD3 monoclonal antibody (M03), clone 6C10

Catalog Number: H00008019-M03

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant BRD3.

Clone Name: 6C10

Immunogen: BRD3 (AAH32124, 418 a.a. ~ 556 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

```
EPVEAPALPAPAAPMVSKGAESSRSSEESSSDSGSSD  
SEEERATRLAELQEQLKAVHEQLAALSQAPVKNPKKK  
KEKKEKEKKKKDKKEKEKEKHKVKAEKKKAKVAPPAK  
QAQKKKAPAKKANSTTTAGRDHFLTCGV
```

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, WB-Ce, WB-Re, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 8019

Gene Symbol: BRD3

Gene Alias: FLJ23227, FLJ41328, KIAA0043, ORFX, RING3L

Gene Summary: This gene was identified based on its homology to the gene encoding the RING3 protein, a

serine/threonine kinase. The gene localizes to 9q34, a region which contains several major histocompatibility complex (MHC) genes. The function of the encoded protein is not known. [provided by RefSeq]