

Datasheet

MKKNK1 monoclonal antibody (M08), clone 2H8

Catalog Number: H00008569-M08

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full length recombinant MKKNK1.

Clone Name: 2H8

Immunogen: MKKNK1 (AAH02755.1, 1 a.a. ~ 465 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MVSSQKLEKPIEMGSSEPLPIADGDRRRKKKRRRGRAT
DSLPGKFEDMYKLTSELLGEGAYAKVQGAVSLQNGKE
YAVKIIKQAGHSRSRVFREVETLYQCQRNKNILELIEF
FEDDTRFYLVFEKLQGGSI LAHIQKQKHFNEREASRVV
RGVAAALDFLHTKDKVSLCHLGWSAMAPSGLTAAPTS
LGSSDPPTSASQVAGTTGIAHRDLKPENILCESPEKVS
PVKICDFDLGSGMKLNNCTPITPELTTPCGSAEYMA
PEVVEVFTDQATFYDKRCDLWSLGVVLYIMLSGYPPF
VGHCGADCGWDRGEVCRVCQNKL FESI QEGKYEFDP
KDWAHISSEAKDLISKLLVRDAKQRLSAAQVLQHPWV
QGQAPEKGLPTPQVLQRNSSTMDLTLFAAEIALNRLL
SQHEENELAEPEALADGLCSMKLSPPCKSRLARRRA
LAQAGRGEDRSPPTAL

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, S-ELISA, WB-Ce, WB-Re, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to

avoid repeated freezing and thawing.

Entrez GeneID: 8569

Gene Symbol: MKKNK1

Gene Alias: MNK1