

Datasheet

NOL3 monoclonal antibody (M02), clone 3A9

Catalog Number: H00008996-M02

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant NOL3.

Clone Name: 3A9

Immunogen: NOL3 (NP_003937.1, 1 a.a. ~ 208 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MGNAQERPSETIDRERKRLVETLQADSGLLLDALLAR
GVLTGPEYEALDALPDAERRVRRLLLLVQKGGEAACQ
ELLRCAQRTAGAPDPAWDWQHVGPGYRDRSYDPPC
PGHWTPEAPGSGTTCPLPRASDPDEAGGPEGSEAV
QSGTPEEPEPELEAEASKEAEPEPEPEPELEPEAEAE
PEPELEPEPDPEPEPDFEERDESEDS

Host: Mouse

Reactivity: Human

Applications: ELISA, S-ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 8996

Gene Symbol: NOL3

Gene Alias: ARC, CARD2, MYC, MYP, NOP, NOP30