

## Datasheet

### **NOL3 monoclonal antibody (M03A), clone 6F5**

**Catalog Number:** H00008996-M03A

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a full-length recombinant NOL3.

**Clone Name:** 6F5

**Immunogen:** NOL3 (NP\_003937.1, 1 a.a. ~ 208 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

MGNAQERPSETIDRERKRLVETLQADSGLLLDALLAR  
GVLTGPEYEALDALPDAERRVRRLLLLVQKGGEAACQ  
ELLRCAQRTAGAPDPAWDWQHVGPGYRDRSYDPPC  
PGHWTPEAPGSGTTCPLPRASDPDEAGGPEGSEAV  
QSGTPEEPEPELEAEASKEAEPEPEPEPELEPEAEAE  
PEPELEPEPDPEPEPDFEERDESEDS

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, WB-Re, WB-Tr  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at  
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG2a Kappa

**Storage Buffer:** In ascites fluid

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 8996

**Gene Symbol:** NOL3

**Gene Alias:** ARC, CARD2, MYC, MYP, NOP, NOP30