

Datasheet

FAIM3 purified MaxPab mouse polyclonal antibody (B02P)

Catalog Number: H00009214-B02P

Regulatory Status: For research use only (RUO)

Product Description: Mouse polyclonal antibody raised against a full-length human FAIM3 protein.

Immunogen: FAIM3 (NP_005440.1, 1 a.a. ~ 390 a.a) full-length human protein.

Sequence:

MDFWLWPLYFLPVSGALRILPEVKVEGELGGSVTIKCP
LPEMHVRIYLCREMAGSGTCGTVVSTTNFIKAEYKGR
VTLKQYPRKNLFLVEVTQLTESDSGVYACGAGMNTDR
GKTQKVTLNVHSEYEPSWEEQPMPETPKWFHLPYLF
QMPAYASSKFVTRVTTTPAQRGKVPVHHSSPTTQIT
HRPRVSRASSVAGDKPRTFLPSTTASKISALEGLLKPO
TPSYNHHTLRHRQALDYGSQSGREGQGFHILPTILG
LFLALLGLVVKRAVERRKALSRRARRLAVRMRALESS
QRPRGSPRPRSQNNIYSACPRRARGADAAGTGEAPV
PGPGAPLPPAPLQVSESPWLHAPSLKTSCEYVSLYHQ
PAAMMEDSDSDDYINVPA

Host: Mouse

Reactivity: Human

Applications: IHC-P, WB-Tr

(See our web site product page for detailed applications information)

Protocols: See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 9214

Gene Symbol: FAIM3

Gene Alias: TOSO