

## Datasheet

### STK17B monoclonal antibody (M01), clone 2C3

**Catalog Number:** H00009262-M01

**Regulatory Status:** For research use only (RUO)

**Product Description:** Mouse monoclonal antibody raised against a partial recombinant STK17B.

**Clone Name:** 2C3

**Immunogen:** STK17B (AAH16040, 272 a.a. ~ 372 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

**Sequence:**

LLVKNPEKRPTAEICLSHSWLQQWDFENLFHPEETS  
SSSQTDHSDVRSSDKTSKSSCNGTCGDREDKENIPE  
DSSMVKRFRFDDSLPNPHELVSDLLC

**Host:** Mouse

**Reactivity:** Human

**Applications:** ELISA, IF, IHC-P, S-ELISA, WB-Ce, WB-Re  
(See our web site product page for detailed applications information)

**Protocols:** See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Isotype:** IgG2b Kappa

**Storage Buffer:** In 1x PBS, pH 7.4

**Storage Instruction:** Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 9262

**Gene Symbol:** STK17B

**Gene Alias:** DRAK2