

Datasheet

TAOK2 monoclonal antibody (M03), clone 2E2

Catalog Number: H00009344-M03

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant TAOK2.

Clone Name: 2E2

Immunogen: TAOK2 (AAH51798, 831 a.a. ~ 930 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

ELQCRQYKRKMLLARHSLDQDLLREDLNKKQTQKDLE
CALLLRQHEATRELELRQLQAVQRTRAELTRLQHQTE
LGNQLEYNKRREQELRQKHAAQVRQQ

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 9344

Gene Symbol: TAOK2

Gene Alias: KIAA0881, MAP3K17, PSK, PSK1, TAO1, TAO2