

Datasheet

MTSS1 monoclonal antibody (M01A), clone 2G9

Catalog Number: H00009788-M01A

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant MTSS1.

Clone Name: 2G9

Immunogen: MTSS1 (AAH23998, 346 a.a. ~ 445 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

GVLPPAPDGPPEERGEHSPESPSVGEGPQGVTSMPS
MWSGQASVNPPLPGPKPSIPEEHRQAIPSEEAEDQER
EPPSATVSPGQIPESDPADLSPRDTPQ

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Ti
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 Kappa

Storage Buffer: In ascites fluid

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 9788

Gene Symbol: MTSS1

Gene Alias: DKFZp781P2223, FLJ44694, KIAA0429, MIM, MIMA, MIMB