

Datasheet

SH2D3C monoclonal antibody (M01), clone 3B2

Catalog Number: H00010044-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant SH2D3C.

Clone Name: 3B2

Immunogen: SH2D3C (NP_005480, 1 a.a. ~ 110 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MTAVGRRCPALGSRGAAGEPEAGSDYVKFSKEYILD
SSPEKLHKELEEELKLSSTDLRSHAWYHGRIPREVSET
LVQRNGDFLIRDLSLGLDYVLCRWRNQLHFKI

Host: Mouse

Reactivity: Human

Applications: ELISA, IHC-P, S-ELISA, WB-Ce, WB-Re
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG1 kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10044

Gene Symbol: SH2D3C

Gene Alias: CHAT, FLJ39664, NSP3, PRO34088