

Datasheet

RBM5 monoclonal antibody (M01), clone 2B6

Catalog Number: H00010181-M01

Regulatory Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a partial recombinant RBM5.

Clone Name: 2B6

Immunogen: RBM5 (AAH02957, 75 a.a. ~ 184 a.a) partial recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

GYHSDGDYGEHDYRHDISDERESKTIMLRGLPITITES
DIREMMESFEGPQPADVRLMKRKTGVSRGFVFEFY
HLQDATSWMEANQKKLVIQGKHIAMHYSNPRPKFED

Host: Mouse

Reactivity: Human

Applications: ELISA, IF, IHC-P, RNAi-Ab, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at <http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Isotype: IgG2a Kappa

Storage Buffer: In 1x PBS, pH 7.4

Storage Instruction: Store at -20°C or lower. Aliquot to avoid repeated freezing and thawing.

Entrez GeneID: 10181

Gene Symbol: RBM5

Gene Alias: FLJ39876, G15, H37, LUCA15, RMB5